



# ***STIC Search Report***

**EIC 1700**

**STIC Database Tracking Number: 130741**

**TO: Alex Noguera  
Location: REM 8C65  
Art Unit : 1753  
September 7, 2004**

**Case Serial Number: 10/019220**

**From: Kathleen Fuller  
Location: EIC 1700  
REMSSEN 4B28  
Phone: 571/272-2505  
Kathleen.Fuller@uspto.gov**

## **Search Notes**

130741

U.S. DEPARTMENT OF COMMERCE  
Patent and Trademark Office

## SEARCH REQUEST FORM

Requestor's

Name:

Alex Noguera

Serial

Number: 10/019,220

Date:

8/25/04

Phone:

571 272-1343

Art Unit:

1753

## Search Topic:

Please write a detailed statement of search topic. Describe specifically as possible the subject matter to be searched. Define any terms that may have a special meaning. Give examples or relevant citations, authors keywords, etc., if known. For sequences, please attach a copy of the sequence. You may include a copy of the broadest and/or most relevant claim(s).

## STAFF USE ONLY

Date completed:

Searcher:

K. Fuller

Terminal time:

40

Elapsed time:

CPU time:

Total time:

70

Number of Searches:

1

Number of Databases:

Search Site

1700

☒ SHC☐ CM-1☐ Pre-S

Type of Search

☐ N.A. Sequence☐ A.A. Sequence☒ Structure☐ Bibliographic

Vendors

☐ IG Suite☒ STN☐ Dialog☐ APS☐ Geninfo☐ SDC☐ DARC/Questel☐ Other

130741  
SEARCH REQUEST FORM

Requestor's Name: Alex Noguera Serial Number: 10/019,220  
Date: 8/25/04 Phone: 571 272-1343 Art Unit: 1753  
8C65

## Search Topic:


Please write a detailed statement of search topic. Describe specifically as possible the subject matter to be searched. Define any terms that may have a special meaning. Give examples or relevant citations, authors keywords, etc., if known. For sequences, please attach a copy of the sequence. You may include a copy of the broadest and/or most relevant claim(s).

22. New) An amperometric sensor suitable for determining the concentration of hydrogen peroxide in a sample, said sensor comprising a ferricyanide compound which, in reduced form, functions as a mediator selective for hydrogen peroxide.

AUG 25


Pat. &amp; T.M. Office

25. (New) A sensor according to claim 22, wherein the ferricyanide compound is of general formula:



in which the groups X are the same or different and at least one X is a non-metallic ion.

40. A ferricyanide compound of formula:



in which the groups X are the same or different and each is a quaternary ammonium ion, at least one of the quaternary ammonium ions having (a) four identical alkyl groups of 5 to 11 carbon atoms other than heptyl or (b) three methyl groups and an alkyl group of 6 to 20 carbon atoms other than hexadecyl.

## STAFF USE ONLY

Date completed: \_\_\_\_\_  
Searcher: K. Fuller  
Terminal time: 40  
Elapsed time: \_\_\_\_\_  
CPU time: \_\_\_\_\_  
Total time: 70  
Number of Searches: 1  
Number of Databases: \_\_\_\_\_

Search Site 1700  
✓ STIG  
\_\_\_\_\_ CM-1  
\_\_\_\_\_ Pre-S

Type of Search  
\_\_\_\_\_ N.A. Sequence  
\_\_\_\_\_ A.A. Sequence  
1 Structure  
\_\_\_\_\_ Bibliographic

Vendors  
\_\_\_\_\_ IG Suite  
✓ STN  
\_\_\_\_\_ Dialog  
\_\_\_\_\_ APS  
\_\_\_\_\_ Geninfo  
\_\_\_\_\_ SDC  
\_\_\_\_\_ DARC/Questel  
\_\_\_\_\_ Other



# STIC Search Results Feedback Form

**EIC17000**

Questions about the scope or the results of the search? Contact *the EIC searcher* or contact:

Kathleen Fuller, EIC 1700 Team Leader  
571/272-2505 REMSEN 4B28

## Voluntary Results Feedback Form

- I am an examiner in Workgroup:  Example: 1713  
➤ Relevant prior art **found**, search results used as follows:

- ☐ 102 rejection
- ☐ 103 rejection
- ☐ Cited as being of interest.
- ☐ Helped examiner better understand the invention.
- ☐ Helped examiner better understand the state of the art in their technology.

Types of relevant prior art found:


- ☐ Foreign Patent(s)
- ☐ Non-Patent Literature  
(journal articles, conference proceedings, new product announcements etc.)

- Relevant prior art **not found**:

- ☐ Results verified the lack of relevant prior art (helped determine patentability).
- ☐ Results were not useful in determining patentability or understanding the invention.

Comments:

Drop off or send completed forms to EIC1700 REMSEN 4B28



=> FILE REG

FILE 'REGISTRY' ENTERED AT 16:47:53 ON 07 SEP 2004  
USE IS SUBJECT TO THE TERMS OF YOUR STN CUSTOMER AGREEMENT.  
PLEASE SEE "HELP USAGETERMS" FOR DETAILS.  
COPYRIGHT (C) 2004 American Chemical Society (ACS)

Property values tagged with IC are from the ZIC/VINITI data file  
provided by InfoChem.

STRUCTURE FILE UPDATES: 6 SEP 2004 HIGHEST RN 740796-45-6  
DICTIONARY FILE UPDATES: 6 SEP 2004 HIGHEST RN 740796-45-6

TSCA INFORMATION NOW CURRENT THROUGH MAY 21, 2004

Please note that search-term pricing does apply when  
conducting SmartSELECT searches.

Crossover limits have been increased. See HELP CROSSOVER for details.

Experimental and calculated property data are now available. For more  
information enter HELP PROP at an arrow prompt in the file or refer  
to the file summary sheet on the web at:  
<http://www.cas.org/ONLINE/DBSS/registryss.html>

=> FILE HCAPLUS

FILE 'HCAPLUS' ENTERED AT 16:47:58 ON 07 SEP 2004  
USE IS SUBJECT TO THE TERMS OF YOUR STN CUSTOMER AGREEMENT.  
PLEASE SEE "HELP USAGETERMS" FOR DETAILS.  
COPYRIGHT (C) 2004 AMERICAN CHEMICAL SOCIETY (ACS)

Copyright of the articles to which records in this database refer is  
held by the publishers listed in the PUBLISHER (PB) field (available  
for records published or updated in Chemical Abstracts after December  
26, 1996), unless otherwise indicated in the original publications.  
The CA Lexicon is the copyrighted intellectual property of the  
the American Chemical Society and is provided to assist you in searching  
databases on STN. Any dissemination, distribution, copying, or storing  
of this information, without the prior written consent of CAS, is  
strictly prohibited.

FILE COVERS 1907 - 7 Sep 2004 VOL 141 ISS 11  
FILE LAST UPDATED: 6 Sep 2004 (20040906/ED)

This file contains CAS Registry Numbers for easy and accurate  
substance identification.

=> D QUE L39

L29 STR

4  
Ak  
2  
Ak~N~Ak  
1 3  
Ak  
5

NODE ATTRIBUTES:


DEFAULT MLEVEL IS ATOM  
DEFAULT ECLEVEL IS LIMITED  
ECOUNT IS M5 C AT 3

GRAPH ATTRIBUTES:

RING(S) ARE ISOLATED OR EMBEDDED  
NUMBER OF NODES IS 5

STEREO ATTRIBUTES: NONE

L32 STR



*25 structures from  
the query*

NODE ATTRIBUTES:

DEFAULT MLEVEL IS ATOM  
DEFAULT ECLEVEL IS LIMITED

GRAPH ATTRIBUTES:

RING(S) ARE ISOLATED OR EMBEDDED  
NUMBER OF NODES IS 7

STEREO ATTRIBUTES: NONE

L34 25 SEA FILE=REGISTRY SSS FUL L32 AND L29  
L36 1 SEA FILE=REGISTRY ABB=ON HYDROGEN PEROXIDE/CN  
L37 13 SEA FILE=HCAPLUS ABB=ON L34  
L38 82929 SEA FILE=HCAPLUS ABB=ON L36  
L39 1 SEA FILE=HCAPLUS ABB=ON L37 AND L38

*13 CA reference*

*only 1 CA references  
with H<sub>2</sub>O<sub>2</sub>*

=> D L39 BIB ABS IND HITSTR

L39 ANSWER 1 OF 1 HCAPLUS COPYRIGHT 2004 ACS on STN

AN 2001:12656 HCAPLUS

DN 134:65598

TI Amperometric sensor for hydrogen peroxide and glucose determination

IN Lau, Kim King Tong; Slater, Jonathan Mark

PA Drew Scientific Limited, UK

SO PCT Int. Appl., 21 pp.

CODEN: PIXXD2

DT Patent

LA English

FAN.CNT 1

*applicant*

|    | PATENT NO.   | KIND | DATE     | APPLICATION NO. | DATE     |
|----|--|------|----------|-----------------|----------|
| PI | WO 2001000865  | A2   | 20010104 | WO 2000-GB2504  | 20000629 |
|    | WO 2001000865  | A3   | 20010913 |                 |          |
|    | W: CA, US  |      |          |                 |          |
|    | RW: AT, BE, CH, CY, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE |      |          |                 |          |
|    | EP 1194585   | A2   | 20020410 | EP 2000-940660  | 20000629 |

KATHLEEN FULLER EIC 1700 REMSEN 4B28 571/272-2505

R: AT, BE, CH, DE, DK, ES, FR, GB, GR, IT, LI, LU, NL, SE, MC, PT,  
IE, FI

PRAI GB 1999-15181 A 19990629  
WO 2000-GB2504 W 20000629

OS MARPAT 134:65598

AB An amperometric sensor suitable for determining the concentration of hydrogen peroxide

in a sample, said sensor comprising a ferricyanide compound which, in reduced form, functions as a mediator specific to hydrogen peroxide.

IC ICM C12Q001-00

CC 79-2 (Inorganic Analytical Chemistry)

ST hydrogen peroxide detn amperometric sensor; glucose sensor

IT Polyamides, analysis

RL: ARU (Analytical role, unclassified); DEV (Device component use); ANST (Analytical study); USES (Uses)

(acrylic; glucose determination by amperometric sensor with ferricyanide compound

as mediator and glucose oxidase)

IT Sensors

(amperometric; hydrogen peroxide determination by amperometric sensor with ferricyanide compound as mediator and glucose oxidase)

IT Glucose sensors

(glucose determination by amperometric sensor with ferricyanide compound as mediator and glucose oxidase)

IT Phosphonium compounds

Quaternary ammonium compounds, analysis

RL: ARU (Analytical role, unclassified); DEV (Device component use); ANST (Analytical study); USES (Uses)

(glucose determination by amperometric sensor with ferricyanide compound as mediator and glucose oxidase)

IT Acrylic polymers, analysis

RL: ARU (Analytical role, unclassified); DEV (Device component use); ANST (Analytical study); USES (Uses)

(polyamide-; glucose determination by amperometric sensor with ferricyanide compound as mediator and glucose oxidase)

IT 50-99-7, Glucose, analysis

RL: ANT (Analyte); ANST (Analytical study)

(glucose determination by amperometric sensor with ferricyanide compound as mediator and glucose oxidase)

IT 79-06-1D, Acrylamide, polymeric derivs. quaterized ferricyanide salts

110-86-1D, Pyridine, polymeric derivs. quaterized ferricyanide salts,

analysis 55066-68-7 58375-66-9 313511-66-9

313511-73-8 313511-84-1 313511-88-5

313511-91-0 313511-94-3 313511-97-6

RL: ARU (Analytical role, unclassified); DEV (Device component use); ANST (Analytical study); USES (Uses)

(glucose determination by amperometric sensor with ferricyanide compound as mediator and glucose oxidase)

IT 313511-61-4P

RL: ARU (Analytical role, unclassified); DEV (Device component use); PNU (Preparation, unclassified); ANST (Analytical study); PREP (Preparation); USES (Uses)

(glucose determination by amperometric sensor with ferricyanide compound as mediator and glucose oxidase)

IT 4328-13-6, Tetrahexylammonium bromide 13746-66-2, Potassium ferricyanide

RL: RCT (Reactant); RACT (Reactant or reagent)

(glucose determination by amperometric sensor with ferricyanide compound as mediator and glucose oxidase)

IT 7722-84-1, Hydrogen peroxide, analysis  
 RL: ANT (Analyte); ANST (Analytical study)  
 (hydrogen peroxide determination by amperometric sensor with ferricyanide compound as mediator and glucose oxidase)

IT 9001-37-0, Glucose oxidase  
 RL: ARG (Analytical reagent use); DEV (Device component use); ANST (Analytical study); USES (Uses)  
 (hydrogen peroxide determination by amperometric sensor with ferricyanide compound as mediator and glucose oxidase)


IT 55066-68-7 313511-66-9 313511-73-8  
 313511-84-1 313511-88-5 313511-91-0  
 313511-94-3 313511-97-6  
 RL: ARU (Analytical role, unclassified); DEV (Device component use); ANST (Analytical study); USES (Uses)  
 (glucose determination by amperometric sensor with ferricyanide compound as mediator and glucose oxidase)

RN 55066-68-7 HCAPLUS

CN 1-Hexadecanaminium, N,N,N-trimethyl-, (OC-6-11)-hexakis(cyano-κC)ferrate(3-) (3:1) (9CI) (CA INDEX NAME)

CM 1

CRN 13408-62-3  
 CMF C6 Fe N6  
 CCI CCS



CM 2

CRN 6899-10-1  
 CMF C19 H42 N

Me<sub>3</sub><sup>+</sup>N- (CH<sub>2</sub>)<sub>15</sub>-Me


RN 313511-66-9 HCAPLUS

CN 1-Decanaminium, N,N,N-tris(decyl)-, (OC-6-11)-hexakis(cyano-κC)ferrate(3-) (3:1) (9CI) (CA INDEX NAME)

CM 1

CRN 48078-03-1  
 CMF C40 H84 N






CM 2

CRN 13408-62-3

CMF C6 Fe N6

CCI CCS



RN 313511-73-8 HCAPLUS


CN 1-Tetradecanaminium, N,N,N-trimethyl-, (OC-6-11)-hexakis(cyano-κC)ferrate(3-) (3:1) (9CI) (CA INDEX NAME)

CM 1

CRN 13408-62-3

CMF C6 Fe N6


CCI CCS



CM 2

CRN 10182-92-0

CMF C17 H38 N




RN 313511-84-1 HCAPLUS

CN 1-Hexanaminium, N,N,N-trimethyl-, (OC-6-11)-hexakis(cyano-  
κC)ferrate(3-) (3:1) (9CI) (CA INDEX NAME)

CM 1

CRN 16208-27-8

CMF C9 H22 N




CM 2

CRN 13408-62-3

CMF C6 Fe N6

CCI CCS




RN 313511-88-5 HCAPLUS

CN 1-Heptanaminium, N,N,N-trimethyl-, (OC-6-11)-hexakis(cyano-  
κC)ferrate(3-) (3:1) (9CI) (CA INDEX NAME)

CM 1

CRN 17077-60-0

CMF C10 H24 N




CM 2

CRN 13408-62-3

CMF C6 Fe N6

CCI CCS



RN 313511-91-0 HCAPLUS  
 CN 1-Octanamminium, N,N,N-trimethyl-, (OC-6-11)-hexakis(cyano-κC)ferrate(3-) (3:1) (9CI) (CA INDEX NAME)


CM 1

CRN 15461-38-8  
 CMF C11 H26 N

Me⁻ (CH₂)₇-N⁺Me₃

CM 2

CRN 13408-62-3  
 CMF C6 Fe N6  
 CCI CCS



RN 313511-94-3 HCAPLUS  
 CN 1-Nonanamminium, N,N,N-trimethyl-, (OC-6-11)-hexakis(cyano-κC)ferrate(3-) (3:1) (9CI) (CA INDEX NAME)

CM 1

CRN 35819-23-9  
 CMF C12 H28 N


Me⁻ (CH₂)₈-N⁺Me₃

CM 2

CRN 13408-62-3

CMF C6 Fe N6

CCI CCS




RN 313511-97-6 HCAPLUS

CN 1-Decanaminium, N,N,N-trimethyl-, (OC-6-11)-hexakis(cyano-κC)ferrate(3-) (3:1) (9CI) (CA INDEX NAME)

CM 1

CRN 15053-09-5

CMF C13 H30 N




CM 2

CRN 13408-62-3

CMF C6 Fe N6

CCI CCS



IT 313511-61-4P

RL: ARU (Analytical role, unclassified); DEV (Device component use); PNU (Preparation, unclassified); ANST (Analytical study); PREP (Preparation); USES (Uses)

(glucose determination by amperometric sensor with ferricyanide compound as mediator and glucose oxidase)


RN 313511-61-4 HCAPLUS

CN 1-Hexanaminium, N,N,N-trihexyl-, (OC-6-11)-hexakis(cyano-  
κC)ferrate(3-) (3:1) (9CI) (CA INDEX NAME)

CM 1

CRN 20256-54-6

CMF C24 H52 N




CM 2

CRN 13408-62-3

CMF C6 Fe N6

CCI CCS



IT 7722-84-1, Hydrogen peroxide, analysis

RL: ANT (Analyte); ANST (Analytical study)

(hydrogen peroxide determination by amperometric sensor with ferricyanide compound as mediator and glucose oxidase)

RN 7722-84-1 HCAPLUS

CN Hydrogen peroxide (H2O2) (9CI) (CA INDEX NAME)

HO-OH

=> => D QUE

L12 4 SEA FILE=REGISTRY ABB=ON ("C6FEN6.(C2H3CL)X.2ZN"/MF OR  
"C6FEN6.(C4H5CL)X.2ZN"/MF OR C6FEN6.1/2C2H8N2.3/2CO.H/MF OR  
C6FEN6.1/2C4H3BRS.3C3H9SN.H2O/MF)  
L13 4 SEA FILE=REGISTRY ABB=ON (C6FEN6.1/2C4H3NO2S.3C3H9SN.3H2O/MF  
OR C6FEN6.1/2C4H4S.3C3H9SN.H2O/MF OR C6FEN6.1/2C6FEN6.2H5O2.3H/

KATHLEEN FULLER EIC 1700 REMSEN 4B28 571/272-2505

MF OR C6FEN6.1/2C6FEN6.3/2FE.1/2K/MF)

L14 2 SEA FILE=REGISTRY ABB=ON (C6FEN6.1/2H12N4PT/MF OR C6FEN6.1/2H12N4PT.2NA/MF)

L15 1 SEA FILE=REGISTRY ABB=ON C6FEN6.1/3CS.4/3NI/MF

L16 1 SEA FILE=REGISTRY ABB=ON C6FEN6.14/9H4N.11/9MN/MF

L17 1 SEA FILE=REGISTRY ABB=ON C6FEN6.2/3C4H12N.5/3O2U/MF

L18 1 SEA FILE=REGISTRY ABB=ON C6FEN6.2/3H4N.4/3K.2/3TM/MF

L19 2 SEA FILE=REGISTRY ABB=ON (C6FEN6.2/5H16MO6O22.4/5H/MF OR C6FEN6.22/15BI.2/5NO3/MF)

L20 2 SEA FILE=REGISTRY ABB=ON (C6FEN6.2AG.2H4N.XH2O/MF OR C6FEN6.2AG.H4N/MF)

L21 5 SEA FILE=REGISTRY ABB=ON (C6FEN6.2BI.2NO3/MF OR C6FEN6.2C2H10CUN2O/MF OR C6FEN6.2C2H10CUN2O.4H2O/MF OR C6FEN6.2C2H2N2NI2S2/MF)

L22 8 SEA FILE=REGISTRY ABB=ON (C6FEN6.2C2H8N2.4H/MF OR C6FEN6.2C2H8N2.XH2O.4H/MF OR C6FEN6.2C3H10CUN2.3H2O.HO/MF OR C6FEN6.2C3H10CUN2.CL.5H2O.K/MF OR C6FEN6.2C3H10CUN2.CL.K/MF OR C6FEN6.2C3H10CUN2.HO/MF OR C6FEN6.2C3H10N2.4H/MF OR C6FEN6.2C3H11OPB.2C3H9PB/MF)


L23 8 SEA FILE=REGISTRY ABB=ON (C6FEN6.2C3H9S.2H/MF OR C6FEN6.2C3H9SN.C2H6SN/MF OR C6FEN6.2C3H9SN.C2H6SN.3H2O/MF OR C6FEN6.2C4H12CUN2/MF OR C6FEN6.2C4H12N.2H/MF OR C6FEN6.2C4H12N.2H2O.H/MF OR C6FEN6.2C4H12N.2H2O.LI/MF OR C6FEN6.2C4H12N.3/2H2O.2H/MF)

L24 8 SEA FILE=REGISTRY ABB=ON (C6FEN6.2C4H12N.H/MF OR C6FEN6.2C4H12N.H2O.NA/MF OR C6FEN6.2C4H12N.K/MF OR C6FEN6.2C4H12N.LI/MF OR C6FEN6.2C4H12N.NA/MF OR C6FEN6.2C4H12N.RB/MF OR C6FEN6.2C4H12N.TL/MF OR C6FEN6.2C4H13N3.4H/MF)

L25 8 SEA FILE=REGISTRY ABB=ON (C6FEN6.2C4H14CUN2O2/MF OR C6FEN6.2C4H16CUN4/MF OR C6FEN6.2C4H16CUN4.5H2O/MF OR C6FEN6.2C4H16N4NI.3H2O.NO3/MF OR C6FEN6.2C4H16N4NI.BF4/MF OR C6FEN6.2C4H16N4NI.CLO4/MF OR C6FEN6.2C4H16N4NI.F6P/MF OR C6FEN6.2C4H16N4NI.NO3/MF)

L26 55 SEA FILE=REGISTRY ABB=ON (L12 OR L13 OR L14 OR L15 OR L16 OR L17 OR L18 OR L19 OR L20 OR L21 OR L22 OR L23 OR L24 OR L25)

L29 STR



NODE ATTRIBUTES:

DEFAULT MLEVEL IS ATOM

DEFAULT ECLEVEL IS LIMITED

ECOUNT IS M5 C AT 3


GRAPH ATTRIBUTES:

RING(S) ARE ISOLATED OR EMBEDDED

NUMBER OF NODES IS 5

STEREO ATTRIBUTES: NONE

L32 STR



NODE ATTRIBUTES:  
 DEFAULT MLEVEL IS ATOM  
 DEFAULT ECLEVEL IS LIMITED

GRAPH ATTRIBUTES:  
 RING(S) ARE ISOLATED OR EMBEDDED  
 NUMBER OF NODES IS 7

STEREO ATTRIBUTES: NONE

L34 25 SEA FILE=REGISTRY SSS FUL L32 AND L29  
 L36 1 SEA FILE=REGISTRY ABB=ON HYDROGEN PEROXIDE/CN  
 L37 13 SEA FILE=HCAPLUS ABB=ON L34  
 L38 82929 SEA FILE=HCAPLUS ABB=ON L36  
 L39 1 SEA FILE=HCAPLUS ABB=ON L37 AND L38  
 L40 387 SEA FILE=HCAPLUS ABB=ON L38 AND ?FERRICYANIDE?  
 L41 17 SEA FILE=HCAPLUS ABB=ON L40 AND SENSOR?  
 L42 48 SEA FILE=HCAPLUS ABB=ON L26  
 L43 0 SEA FILE=HCAPLUS ABB=ON L38 AND L42  
 L44 17 SEA FILE=HCAPLUS ABB=ON L41 OR L43  
 L45 16 SEA FILE=HCAPLUS ABB=ON L44 NOT L39

*76 more CA  
 references with*

=> D L45 BIB ABS IND HITSTR 1-16

L45 ANSWER 1 OF 16 HCAPLUS COPYRIGHT 2004 ACS on STN  
 AN 2004:375567 HCAPLUS  
 DN 140:367826  
 TI Printed circuit board electrochemical **sensor**  
 IN Shiu, Tian-Tsai; Jang, Jing-Yu; Wang, Ji-Wen  
 PA Industrial Technology Research Institute, Taiwan  
 SO Taiwan, 4 pp.  
 CODEN: TWXXA5  
 DT Patent  
 LA Chinese  
 FAN.CNT 1

|      | PATENT NO.       | KIND | DATE     | APPLICATION NO.  | DATE     |
|------|------------------|------|----------|------------------|----------|
| PI   | TW 496110        | B    | 20020721 | TW 1998-87117851 | 19981028 |
| PRAI | TW 1998-87117851 |      | 19981028 |                  |          |

AB This invention provides a novel manufacture process of printed circuit board (PCB) electrochem. **sensor**. A precious metal material of required thickness is plated onto an outer surface of a standard PCB to fully cover the conductive circuit substrate of the PCB. According to the invention, all the conductive nodes of the electrodes on the same substrate are linked to an electrode during circuit layout. When preparing the electrochem. electrode, the precious metal electrode material of required thickness is then plated onto an outer surface of a standard PCB to fully cover the conductive circuit substrate of the PCB. There is no substrate

exposed on the cross-section surface to achieve the purpose of using the manufacture and material of PCB for electrochem. electrode. The invention also discloses an electrochem. **sensor** obtained through the manufacture process.

IC ICM H05K003-00  
 CC 79-2 (Inorganic Analytical Chemistry)  
 Section cross-reference(s): 9, 80  
 ST printed circuit electrochem **sensor**  
 IT **Sensors**  
     (electrochem.; printed circuit board electrochem. **sensor**)  
 IT Blood analysis  
     (glucose; printed circuit board electrochem. **sensor**)  
 IT Composites  
     Electric circuits  
     Electric conductivity  
     Electrodes  
     Glass substrates  
     Glucose **sensors**  
     Materials  
     Printed circuit boards  
     Surface  
     Thickness  
     (printed circuit board electrochem. **sensor**)  
 IT Enzymes, uses  
     RL: ARG (Analytical reagent use); DEV (Device component use); ANST  
     (Analytical study); USES (Uses)  
     (printed circuit board electrochem. **sensor**)  
 IT Noble metals  
     RL: DEV (Device component use); USES (Uses)  
     (printed circuit board electrochem. **sensor**)  
 IT Ceramics  
     (substrates; printed circuit board electrochem. **sensor**)  
 IT 50-99-7, Glucose, analysis **7722-84-1**, Hydrogen peroxide,  
     analysis 13746-66-2, Potassium **ferricyanide**  
     RL: ANT (Analyte); ANST (Analytical study)  
     (printed circuit board electrochem. **sensor**)  
 IT 7440-02-0, Nickel, uses 7440-05-3, Palladium, uses 7440-06-4,  
     Platinum, uses 7440-22-4, Silver, uses 7440-50-8, Copper, uses  
     7440-57-5, Gold, uses  
     RL: DEV (Device component use); USES (Uses)  
     (printed circuit board electrochem. **sensor**)  
 IT **7722-84-1**, Hydrogen peroxide, analysis  
     RL: ANT (Analyte); ANST (Analytical study)  
     (printed circuit board electrochem. **sensor**)  
 RN 7722-84-1 HCAPLUS  
 CN Hydrogen peroxide (H2O2) (9CI) (CA INDEX NAME)

HO-OH

L45 ANSWER 2 OF 16 HCAPLUS COPYRIGHT 2004 ACS on STN  
 AN 2003:571175 HCAPLUS  
 DN 139:110751  
 TI Method and apparatus for processing electrochemical signals  
 IN Iyengar, Sridhar G.; Haas, Daniel; Bolon, Craig  
 PA Agamatrix, Inc., USA

KATHLEEN FULLER EIC 1700 REMSEN 4B28 571/272-2505



SO PCT Int. Appl., 72 pp.

CODEN: PIXXD2

DT Patent

LA English

FAN.CNT 1

|      | PATENT NO.  | KIND | DATE     | APPLICATION NO. | DATE     |
|------|---|------|----------|-----------------|----------|
| PI   | WO 2003060154   | A2   | 20030724 | WO 2003-US1113  | 20030115 |
|      | WO 2003060154   | A3   | 20040805 |                 |          |
|      | W: AE, AG, AL, AM, AT, AU, AZ, BA, BB, BG, BR, BY, BZ, CA, CH, CN, CO, CR, CU, CZ, DE, DK, DM, DZ, EC, EE, ES, FI, GB, GD, GE, GH, GM, HR, HU, ID, IL, IN, IS, JP, KE, KG, KP, KR, KZ, LC, LK, LR, LS, LT, LU, LV, MA, MD, MG, MK, MN, MW, MX, MZ, NO, NZ, OM, PH, PL, PT, RO, RU, SC, SD, SE, SG, SK, SL, TJ, TM, TN, TR, TT, TZ, UA, UG, US, UZ, VC, VN, YU, ZA, ZM, ZW, AM, AZ, BY, KG, KZ, MD, RU, TJ, TM<br>RW: GH, GM, KE, LS, MW, MZ, SD, SL, SZ, TZ, UG, ZM, ZW, AT, BE, BG, CH, CY, CZ, DE, DK, EE, ES, FI, FR, GB, GR, HU, IE, IT, LU, MC, NL, PT, SE, SI, SK, TR, BF, BJ, CF, CG, CI, CM, GA, GN, GQ, GW, ML, MR, NE, SN, TD, TG |      |          |                 |          |
|      | US 2003178322   | A1   | 20030925 | US 2003-342794  | 20030115 |
| PRAI | US 2002-350175P   | P    | 20020115 |                 |          |

AB Systems and methods are provided herein for improving the selectivity and productivity of **sensors** via digital signal processing techniques. According to one illustrative embodiment, in an electrochem. method for monitoring of a select analyte in a mixed sample with an interfering analyte, an improvement is provided that includes applying a large amplitude potential stimulus waveform to the sample to generate a nonlinear current signal; and resolving a signal contribution from the select analyte in the generated signal by a vector projection method with an analyte vector comprising a plurality of real and imaginary parts of one or more Fourier coeffs. at one or more frequencies of a reference current signal for the select analyte.

IC ICM C12Q001-00

CC 79-2 (Inorganic Analytical Chemistry)

Section cross-reference(s): 72, 80

ST app electrochem signal processing

IT Mathematical methods

(Fourier-transform; analyte determination in mixts. by electrochem. anal.

and

method and apparatus for processing electrochem. signals for improved selectivity)

IT Polarography

(a.c.; analyte determination in mixts. by electrochem. anal. and method and apparatus for processing electrochem. signals for improved selectivity)

IT Cyclic voltammetry

Data processing

Linear-sweep voltammetry

Potentiostats

Square wave voltammetry

(analyte determination in mixts. by electrochem. anal. and method and apparatus for

processing electrochem. signals for improved selectivity)

IT Enzymes, analysis

RL: ANT (Analyte); ANST (Analytical study)

(analyte determination in mixts. by electrochem. anal. and method and apparatus for

processing electrochem. signals for improved selectivity)

IT Electric circuits  
 (galvanostats; analyte determination in mixts. by electrochem. anal. and method  
 and apparatus for processing electrochem. signals for improved selectivity)  
 IT Electric impedance  
 (spectroscopy; analyte determination in mixts. by electrochem. anal. and method  
 and apparatus for processing electrochem. signals for improved selectivity)  
 IT 50-99-7, D-Glucose, analysis 51-61-6, Dopamine, analysis 69-93-2, Uric acid, analysis 102-54-5, Ferrocene 103-90-2, Acetaminophen 7722-84-1, Hydrogen peroxide, analysis 13408-62-3, **Ferricyanide** 13408-63-4, Ferrocyanide  
 RL: ANT (Analyte); ANST (Analytical study)  
 (analyte determination in mixts. by electrochem. anal. and method and apparatus for processing electrochem. signals for improved selectivity)  
 IT 50-81-7, Ascorbic acid, analysis  
 RL: ARU (Analytical role, unclassified); ANST (Analytical study)  
 (analyte determination in mixts. by electrochem. anal. and method and apparatus for processing electrochem. signals for improved selectivity)  
 IT 7722-84-1, Hydrogen peroxide, analysis  
 RL: ANT (Analyte); ANST (Analytical study)  
 (analyte determination in mixts. by electrochem. anal. and method and apparatus for processing electrochem. signals for improved selectivity)  
 RN 7722-84-1 HCAPLUS  
 CN Hydrogen peroxide (H2O2) (9CI) (CA INDEX NAME)

HO-OH

L45 ANSWER 3 OF 16 HCAPLUS COPYRIGHT 2004 ACS on STN  
 AN 2003:300042 HCAPLUS  
 DN 139:116295  
 TI Measurement system of low glucose concentration during the cultivation of yeast cells  
 AU Kishimoto, Tomokazu; Hara, Seiichi; Muraji, Masafumi; Tsujimoto, Hiroaki; Azuma, Masayuki; Ooshima, Hiroshi  
 CS Department of Physical Electronics and Information, Japan  
 SO Memoirs of the Faculty of Engineering, Osaka City University (2002), 43, 19-23  
 CODEN: MFEOAR; ISSN: 0078-6659  
 PB Osaka City University, Faculty of Engineering  
 DT Journal  
 LA English  
 AB A yeast cell changes an active state in accordance with glucose concentration in  
 a culture medium. Below a certain critical glucose concentration under aerobic conditions, the yeast respire. Exceeding its value, the yeast changes an active state to fermentation. The aim of our study is to maintain the state of respiration and fermentation of yeast artificially. And so, a glucose **sensor** was needed to satisfy with respiration condition. In this study, we tried to construct the glucose **sensor** which was to measure glucose concns. in very low region for a long time and to maintain quasi real-time measurement. The **sensor** was constructed using

the phenomena of light emission by luminol, we evaluated the sensitivity, stability and reliability of it. The **sensor** was robust against outer disturbances, and had an influence by flow rate of solution, and dialysis rate. A detailed explanation of aerobic conditions and of reaction principle of the constructed glucose **sensor** will be presented here. And then, some basic characteristics of the glucose **sensor** will be shown here as well.

CC 16-1 (Fermentation and Bioindustrial Chemistry)  
 ST yeast fermn glucose measurement  
 IT Metabolism  
     (Crabtree effect; system to measure low glucose concns. during yeast fermns.)  
 IT Fermentation  
     (aerobic; system to measure low glucose concns. during yeast fermns.)  
 IT Biosensors  
     (enzymic; system to measure low glucose concns. during yeast fermns.)  
 IT Process control  
     (online; system to measure low glucose concns. during yeast fermns.)  
 IT Fermentation  
     Saccharomyces cerevisiae  
     (system to measure low glucose concns. during yeast fermns.)  
 IT 50-99-7, Dextrose, analysis  
     RL: ANT (Analyte); BCP (Biochemical process); ANST (Analytical study);  
     BIOL (Biological study); PROC (Process)  
     (system to measure low glucose concns. during yeast fermns.)  
 IT 9001-37-0, Glucose oxidase  
     RL: ARG (Analytical reagent use); BCP (Biochemical process); CAT (Catalyst use); ANST (Analytical study); BIOL (Biological study); PROC (Process);  
     USES (Uses)  
     (system to measure low glucose concns. during yeast fermns.)  
 IT 521-31-3, Luminol 13746-66-2, Potassium **ferricyanide**  
     RL: ARG (Analytical reagent use); RCT (Reactant); ANST (Analytical study);  
     RACT (Reactant or reagent); USES (Uses)  
     (system to measure low glucose concns. during yeast fermns.)  
 IT 7722-84-1P, Hydrogen peroxide, preparation  
     RL: BPN (Biosynthetic preparation); RCT (Reactant); BIOL (Biological study); PREP (Preparation); RACT (Reactant or reagent)  
     (system to measure low glucose concns. during yeast fermns.)  
 IT 7722-84-1P, Hydrogen peroxide, preparation  
     RL: BPN (Biosynthetic preparation); RCT (Reactant); BIOL (Biological study); PREP (Preparation); RACT (Reactant or reagent)  
     (system to measure low glucose concns. during yeast fermns.)  
 RN 7722-84-1 HCAPLUS  
 CN Hydrogen peroxide (H2O2) (9CI) (CA INDEX NAME)

HO-OH

RE.CNT 1 THERE ARE 1 CITED REFERENCES AVAILABLE FOR THIS RECORD  
 ALL CITATIONS AVAILABLE IN THE RE FORMAT

L45 ANSWER 4 OF 16 HCAPLUS COPYRIGHT 2004 ACS on STN  
 AN 2003:114720 HCAPLUS  
 DN 140:24956  
 TI Chemiluminescence microfluidic system **sensor** on a chip for  
 determination of glucose in human serum with immobilized reagents  
 AU Lu, Yi; Zhang, Zhujun; Chen, Funan

KATHLEEN FULLER EIC 1700 REMSEN 4B28 571/272-2505

CS Institute of Analytical Science, Department of Chemistry, Southwest Normal University, Chungking, 400715, Peop. Rep. China

SO Talanta (2003), 59(3), 571-576  
CODEN: TLNTA2; ISSN: 0039-9140

PB Elsevier Science B.V.

DT Journal

LA English

AB A chemiluminescence (CL) biosensor on a chip coupled to microfluidic system is described in this paper. The CL biosensor measured 25+45+5 mm in dimension, was readily produced in anal. laboratory Glucose oxidase (GOD) was immobilized onto controlled-pore glass (CPG) via glutaraldehyde activation and packed into a reservoir. The anal. reagents, including luminol and **ferricyanide**, were electrostatically co-immobilized on an anion-exchange resin. The most characteristic of the biosensor was to introduce the air as the carrier flow instead of the common solution carrier for the first. The glucose was sensed by the CL reaction between hydrogen peroxide produced from the enzymic reaction and CL reagents, which were released from the anion-exchange resin. The proposed method has been successfully applied to the determination of glucose in human serum. The linear range of the glucose concentration was 1.1-110 mM and the detection limit was 0.1 mM (3σ).

CC 9-1 (Biochemical Methods)

ST chemiluminescence biosensor glucose detn blood serum

IT Blood analysis  
Blood serum  
Human  
(chemiluminescence microfluidic biosensor on a chip for determination of glucose in human serum with immobilized enzyme and reagents)

IT Biosensors  
(enzymic, chemiluminescent; chemiluminescence microfluidic biosensor on a chip for determination of glucose in human serum with immobilized enzyme and reagents)

IT 50-99-7, D-Glucose, analysis  
RL: ANT (Analyte); ANST (Analytical study)  
(chemiluminescence microfluidic biosensor on a chip for determination of glucose in human serum with immobilized enzyme and reagents)

IT 7722-84-1, Hydrogen peroxide, uses  
RL: ARG (Analytical reagent use); ANST (Analytical study); USES (Uses)  
(chemiluminescence microfluidic biosensor on a chip for determination of glucose in human serum with immobilized enzyme and reagents)

IT 521-31-3D, Luminol, immobilized 9001-37-0D, Glucose oxidase, immobilized 13408-62-3D, **Ferricyanide**, immobilized  
RL: ARG (Analytical reagent use); DEV (Device component use); ANST (Analytical study); USES (Uses)  
(chemiluminescence microfluidic biosensor on a chip for determination of glucose in human serum with immobilized enzyme and reagents)

IT 7722-84-1, Hydrogen peroxide, uses  
RL: ARG (Analytical reagent use); ANST (Analytical study); USES (Uses)  
(chemiluminescence microfluidic biosensor on a chip for determination of glucose in human serum with immobilized enzyme and reagents)

RN 7722-84-1 HCAPLUS

CN Hydrogen peroxide (H2O2) (9CI) (CA INDEX NAME)

HO-OH

RE.CNT 28 THERE ARE 28 CITED REFERENCES AVAILABLE FOR THIS RECORD  
ALL CITATIONS AVAILABLE IN THE RE FORMAT

L45 ANSWER 5 OF 16 HCAPLUS COPYRIGHT 2004 ACS on STN  
AN 2002:535543 HCAPLUS  
DN 137:269618  
TI Electrochemical preparation, characterization and application of  
electrodes modified with hybrid hexacyanoferrates of copper and cobalt  
AU Cui, Xingpin; Hong, Li; Lin, Xiangqin  
CS Department of Chemistry, University of Science and Technology of China,  
Hefei, 230026, Peop. Rep. China  
SO Journal of Electroanalytical Chemistry (2002), 526(1-2), 115-124  
CODEN: JECHES  
PB Elsevier Science B.V.  
DT Journal  
LA English  
AB Hybrid Cu-Co hexacyanoferrate (CuCoHCF) films were electrodeposited on a  
Pt electrode or a glassy C electrode by cyclic voltammetry and  
characterized by electrochem., XRD, ICP-AES and XPS. The results  
indicated that CuCoHCF was a substitution-type hybrid hexacyanoferrate.  
With the increase of Cu<sup>2+</sup> content in the deposition solution, the Cu<sup>2+</sup>  
content in the films increased correspondingly, while the lattice constant  
of the films decreased gradually. The CuCoHCF modified Pt electrode  
exhibited stable electrochem. responses in a wide pH range of 4-10 and  
permeability for monovalent cations in the order of K<sup>+</sup>>Li<sup>+</sup>>Na<sup>+</sup>>NH<sub>4</sub><sup>+</sup>, both  
of which are different from those of the resp. single component Cu or Co  
hexacyanoferrates. XPS gave direct evidence that the Fe element existed  
as Fe(III) in oxidized films and was reduced to Fe(II) during x-ray  
scanning. K<sup>+</sup> was incorporated into and excluded from CuCoHCF films to  
maintain elec. neutrality during the reduction and oxidation process, resp.  
The  
CuCoHCF modified glassy C electrode exhibited obvious electrocatalytic  
activity towards both reduction and oxidation of H<sub>2</sub>O<sub>2</sub>. When a cathodic  
catalytic  
current was used, the **sensor** exhibited a linear response in a  
H<sub>2</sub>O<sub>2</sub> concentration range of 2.3 × 10<sup>-3</sup>-8.1 × 10<sup>-7</sup> M with a detection  
limit of 6.6 × 10<sup>-8</sup> M The H<sub>2</sub>O<sub>2</sub> **sensor** showed excellent  
stability and anti-interference ability towards O and other easily  
oxidized compds. due to a low applied potential of 0.02 V, which is a  
great merit for further application in the field of biosensors.  
CC 72-2 (Electrochemistry)  
Section cross-reference(s): 67, 79  
ST electrochem prepn electrode modified hybrid copper cobalt hexacyanoferrate  
IT Auger electron spectra  
Chemically modified electrodes  
Cyclic voltammetry  
Electrodeposition  
X-ray diffraction  
X-ray photoelectron spectra  
(electrochem. preparation, characterization and application of electrodes  
modified with hybrid hexacyanoferrates of copper and cobalt)  
IT **Sensors**  
(electrochem., for cations; electrochem. preparation, characterization and  
application of electrodes modified with hybrid hexacyanoferrates of  
copper and cobalt)  
IT Oxidation catalysts  
Redox reaction

Reduction catalysts

(electrochem.; electrochem. preparation, characterization and application of electrodes modified with hybrid hexacyanoferrates of copper and cobalt)

- IT Permeability  
(to cations; electrochem. preparation, characterization and application of electrodes modified with hybrid hexacyanoferrates of copper and cobalt)
- IT 13601-13-3P 15415-49-3P 41754-48-7P  
RL: CAT (Catalyst use); DEV (Device component use); PNU (Preparation, unclassified); PRP (Properties); PREP (Preparation); USES (Uses)  
(electrochem. preparation, characterization and application of electrodes modified with hybrid hexacyanoferrates of copper and cobalt)
- IT 7722-84-1, Hydrogen peroxide, properties  
RL: ANT (Analyte); PRP (Properties); ANST (Analytical study)  
(electrochem. reduction and oxidation of; electrochem. preparation, characterization and application of electrodes modified with hybrid hexacyanoferrates of copper and cobalt)
- IT 7440-06-4, Platinum, uses  
RL: DEV (Device component use); USES (Uses)  
(electrode; electrochem. preparation, characterization and application of electrodes modified with hybrid hexacyanoferrates of copper and cobalt)
- IT 7440-44-0, Carbon, uses  
RL: DEV (Device component use); USES (Uses)  
(glassy, electrode; electrochem. preparation, characterization and application of electrodes modified with hybrid hexacyanoferrates of copper and cobalt)
- IT 7757-79-1, Potassium nitrate, uses  
RL: NUU (Other use, unclassified); USES (Uses)  
(in cobalt copper hexacyanoferrate-modified electrode preparation; electrochem. preparation, characterization and application of electrodes modified with hybrid hexacyanoferrates of copper and cobalt)
- IT 7758-98-7, Copper sulfate, reactions 10141-05-6, Cobalt nitrate  
13746-66-2, Potassium **ferricyanide**  
RL: RCT (Reactant); RACT (Reactant or reagent)  
(in cobalt copper hexacyanoferrate-modified electrode preparation; electrochem. preparation, characterization and application of electrodes modified with hybrid hexacyanoferrates of copper and cobalt)
- IT 7439-93-2D, Lithium, ions, properties 7440-09-7D, Potassium, ions, properties 7440-23-5D, Sodium, ions, properties 14798-03-9, Ammonium, properties  
RL: ANT (Analyte); PRP (Properties); ANST (Analytical study)  
(sensing of; electrochem. preparation, characterization and application of electrodes modified with hybrid hexacyanoferrates of copper and cobalt)
- IT 7722-84-1, Hydrogen peroxide, properties  
RL: ANT (Analyte); PRP (Properties); ANST (Analytical study)  
(electrochem. reduction and oxidation of; electrochem. preparation, characterization and application of electrodes modified with hybrid hexacyanoferrates of copper and cobalt)
- RN 7722-84-1 HCAPLUS
- CN Hydrogen peroxide (H2O2) (9CI) (CA INDEX NAME)

HO-OH

RE.CNT 59 THERE ARE 59 CITED REFERENCES AVAILABLE FOR THIS RECORD  
ALL CITATIONS AVAILABLE IN THE RE FORMAT

L45 ANSWER 6 OF 16 HCAPLUS COPYRIGHT 2004 ACS on STN

KATHLEEN FULLER EIC 1700 REMSEN 4B28 571/272-2505

AN 2000:109096 HCAPLUS  
 DN 132:331556  
 TI Flexible amperometric transducers for biosensors based on a screen printed three electrode system  
 AU Erlenkotter, A.; Kottbus, M.; Chemnitz, G. -C.  
 CS Department of Inorganic Chemistry, University of Munster, Munster, D-48149, Germany  
 SO Journal of Electroanalytical Chemistry (2000), 481(1), 82-94  
 CODEN: JECHES; ISSN: 0368-1874  
 PB Elsevier Science S.A.  
 DT Journal  
 LA English  
 AB Screen printed three electrode **sensors** comprising a platinum working, a carbon counter and an Ag|AgCl pseudo reference electrode were developed employing polymer thick film inks. The **sensors** were constructed as amperometric transducers for multianalyte biosensors for use in batch, as well as in flow through systems. The characteristics of the **sensors** were determined. The active surface area of the Pt working electrodes was determined using electrochem. and SEM studies. Cyclic voltammograms of the **ferricyanide/ferrocyanide** couple showed that the reaction was quasi-reversible at these electrodes. Although the surface was not ideal for this redox couple, the **sensors** proved to be reproducible and well suited for the determination of hydrogen peroxide and thus for biosensors based on oxidases as biol. active compds. The combination of two pretreatment steps, an addnl. heat curing and an electrochem. preconditioning step, was found to be most helpful to reduce background current and settling time of the **sensors**. Different aspects of the changing surface composition are discussed. The **sensors** with optimized preconditioning showed linear ranges from 10  $\mu\text{M}$  up to at least 500  $\mu\text{M}$  hydrogen peroxide and sensitivities of  $6.97 \pm 0.20 \text{ nA } \mu\text{M}^{-1}$  hydrogen peroxide for uncovered,  $4.01 \pm 0.08 \text{ nA } \mu\text{M}^{-1}$  hydrogen peroxide for PCS/BSA membrane covered and  $0.222 \pm 0.002 \text{ nA } \mu\text{M}^{-1}$  hydrogen peroxide for Nafion® coated platinum working electrodes. Moreover, optimized transducers with immobilized sarcosine oxidase (sensitivity:  $2.30 \pm 0.07 \text{ nA } \mu\text{M}^{-1}$  sarcosine) demonstrated the feasibility of the **sensor** concept, the manufacturing and pretreatment processes for the development of enzyme **sensors**.  
 CC 9-7 (Biochemical Methods)  
 Section cross-reference(s): 6, 7, 72  
 ST amperometric electrode biosensor screen printing enzyme  
 IT Biosensors  
 Electrodes  
 Enzyme electrodes  
 (amperometric; flexible amperometric transducers for biosensors based on a screen printed three electrode system)  
 IT Electrodes  
 Screen printing  
 (flexible amperometric transducers for biosensors based on a screen printed three electrode system)  
 IT 7722-84-1, Hydrogen peroxide, analysis  
 RL: ANT (Analyte); ANST (Analytical study)  
 (flexible amperometric transducers for biosensors based on a screen printed three electrode system)  
 IT 107-97-1, Sarcosine  
 RL: ANT (Analyte); BPR (Biological process); BSU (Biological study, unclassified); ANST (Analytical study); BIOL (Biological study); PROC (Process)

(flexible amperometric transducers for biosensors based on a screen printed three electrode system)

IT 9029-22-5, Sarcosine oxidase 9035-73-8, Oxidase  
 RL: ARU (Analytical role, unclassified); BAC (Biological activity or effector, except adverse); BPR (Biological process); BSU (Biological study, unclassified); ANST (Analytical study); BIOL (Biological study); PROC (Process)

(flexible amperometric transducers for biosensors based on a screen printed three electrode system)

IT 7722-84-1, Hydrogen peroxide, analysis  
 RL: ANT (Analyte); ANST (Analytical study)

(flexible amperometric transducers for biosensors based on a screen printed three electrode system)

RN 7722-84-1 HCAPLUS

CN Hydrogen peroxide (H2O2) (9CI) (CA INDEX NAME)

HO-OH

RE.CNT 49 THERE ARE 49 CITED REFERENCES AVAILABLE FOR THIS RECORD  
 ALL CITATIONS AVAILABLE IN THE RE FORMAT

L45 ANSWER 7 OF 16 HCAPLUS COPYRIGHT 2004 ACS on STN  
 AN 2000:47107 HCAPLUS  
 DN 132:87444  
 TI Chemical sensing techniques employing liquid-core optical fibers  
 IN Fein, Harry; Liu, Su-yi  
 PA World Precision Instruments, Inc., USA  
 SO U.S., 13 pp., Cont.-in-part of U. S. Ser. No. 951,254.  
 CODEN: USXXAM

DT Patent  
 LA English

FAN.CNT 2

|      | PATENT NO.  | KIND | DATE     | APPLICATION NO. | DATE     |
|------|---|------|----------|-----------------|----------|
| PI   | US 6016372  | A    | 20000118 | US 1998-55865   | 19980406 |
|      | US 6011882  | A    | 20000104 | US 1997-951254  | 19971016 |
|      | EP 909946   | A2   | 19990421 | EP 1998-308481  | 19981016 |
|      | EP 909946   | A3   | 19990811 |                 |          |
|      | EP 909946   | B1   | 20040623 |                 |          |
|      | R: AT, BE, CH, DE, DK, ES, FR, GB, GR, IT, LI, LU, NL, SE, MC, PT, IE, SI, LT, LV, FI, RO |      |          |                 |          |
| PRAI | US 1997-951254  | A2   | 19971016 |                 |          |
|      | US 1998-55865   | A    | 19980406 |                 |          |

AB A gas or vapor permeable optical fiber waveguide with a liquid core is employed as a probe for the detection or measurement of a chemical specie of interest by filling the waveguide core region with a reagent liquid which undergoes a change in an optical characteristic thereof when exposed to the chemical specie and then inserting the filled waveguide into an environment in which the chemical specie may be present. The chemical specie, if present, will permeate through the waveguide wall and react with or be absorbed in the core liquid. Sensitivity is enhanced by controlling the pressure differential across the waveguide wall and/or by shaping the waveguide to enlarge the surface area. When the reaction generates light, the devices which detect that light will be shaped and disposed to maximize the collection thereof. The **sensor** is suitable for applications including CO2, O2, CO, H2S, NO2 NH3, O3, H2O2, chlorine,



concentrated acids detection; detection of organic compds. in water;  
 respiratory  
 air anal.; and an in-line monitor for control purposes.

IC ICM G02B006-20  
 NCL 385012000  
 CC 79-6 (Inorganic Analytical Chemistry)  
 Section cross-reference(s): 9, 59, 61, 80

ST liq core optical fiber **sensor**  
 IT Absorption spectroscopy  
 Acid-base indicators  
 (carbon dioxide detection by optical **sensor** with  
 liquid-core-filled optical fiber waveguide containing pH indicator in  
 aqueous carbonate-bicarbonate buffer)

IT Optical gas **sensors**  
 Optical gas **sensors**  
 (fiber-optic; gas **sensor** based on optical properties of  
 liquid-core-filled optical fiber waveguides)

IT Fluorometry  
 Gas analysis  
 Luminescence spectroscopy  
 Optical fibers  
 Optical waveguides  
 Raman spectroscopy  
 Respiratory air  
 (gas **sensor** based on optical properties of liquid-core-filled  
 optical fiber waveguides)

IT Fiber optic **sensors**  
 Fiber optic **sensors**  
 (gas; gas **sensor** based on optical properties of  
 liquid-core-filled optical fiber waveguides)

IT Acids, analysis  
 RL: ANT (Analyte); ANST (Analytical study)  
 (inorg.; concentrated inorg. acids detection by optical **sensor** with  
 liquid-core-filled optical fiber waveguide containing indicator solution)

IT Organic compounds, analysis  
 RL: ANT (Analyte); ANST (Analytical study)  
 (organic compds. detection in water by optical **sensor** with  
 liquid-core-filled optical fiber waveguide containing indicator solution)

IT 7664-41-7, Ammonia, analysis  
 RL: ANT (Analyte); ANST (Analytical study)  
 (ammonia detection by optical **sensor** with liquid-core-filled  
 optical fiber waveguide containing indicator solution)

IT 124-38-9, Carbon dioxide, analysis  
 RL: ANT (Analyte); ANST (Analytical study)  
 (carbon dioxide detection by optical **sensor** with  
 liquid-core-filled optical fiber waveguide containing pH indicator in  
 aqueous carbonate-bicarbonate buffer)

IT 76-59-5, Bromothymol blue 143-74-8, Phenol red  
 RL: ARG (Analytical reagent use); DEV (Device component use); ANST  
 (Analytical study); USES (Uses)  
 (carbon dioxide detection by optical **sensor** with  
 liquid-core-filled optical fiber waveguide containing pH indicator in  
 aqueous carbonate-bicarbonate buffer)

IT 630-08-0, Carbon monoxide, analysis  
 RL: ANT (Analyte); ANST (Analytical study)

(carbon monoxide detection by optical **sensor** with liquid-core-filled optical fiber waveguide containing indicator solution)  
IT 10025-98-6, Potassium tetrachloropalladate(II) 222159-57-1, Silver p-sulfoaminobenzoate  
RL: ARG (Analytical reagent use); DEV (Device component use); ANST (Analytical study); USES (Uses)  
(carbon monoxide detection by optical **sensor** with liquid-core-filled optical fiber waveguide containing indicator solution)  
IT 7782-50-5, Chlorine, analysis  
RL: ANT (Analyte); ANST (Analytical study)  
(chlorine detection by optical **sensor** with liquid-core-filled optical fiber waveguide containing indicator solution)  
IT 119-93-7, o-Tolidine 34314-06-2, Tetramethylbenzidine  
RL: ARG (Analytical reagent use); DEV (Device component use); ANST (Analytical study); USES (Uses)  
(chlorine detection by optical **sensor** with liquid-core-filled optical fiber waveguide containing indicator solution)  
IT 7722-84-1, Hydrogen peroxide, analysis  
RL: ANT (Analyte); ANST (Analytical study)  
(hydrogen peroxide detection by optical **sensor** with liquid-core-filled optical fiber waveguide containing indicator solution)  
IT 521-31-3, Luminol  
RL: ARG (Analytical reagent use); DEV (Device component use); ANST (Analytical study); USES (Uses)  
(hydrogen peroxide detection by optical **sensor** with liquid-core-filled optical fiber waveguide containing indicator solution)  
IT 13746-66-2, Potassium **ferricyanide**  
RL: CAT (Catalyst use); DEV (Device component use); USES (Uses)  
(hydrogen peroxide detection by optical **sensor** with liquid-core-filled optical fiber waveguide containing indicator solution)  
IT 7783-06-4, Hydrogen sulfide (H<sub>2</sub>S), analysis  
RL: ANT (Analyte); ANST (Analytical study)  
(hydrogen sulfide detection by optical **sensor** with liquid-core-filled optical fiber waveguide containing indicator solution)  
IT 14402-89-2, Sodium nitroprusside  
RL: ARG (Analytical reagent use); DEV (Device component use); ANST (Analytical study); USES (Uses)  
(hydrogen sulfide detection by optical **sensor** with liquid-core-filled optical fiber waveguide containing indicator solution)  
IT 10102-44-0, Nitrogen oxide (NO<sub>2</sub>), analysis  
RL: ANT (Analyte); ANST (Analytical study)  
(nitrogen dioxide detection by optical **sensor** with liquid-core-filled optical fiber waveguide containing indicator solution)  
IT 63-74-1, Sulfanilamide 121-57-3, Sulfanilic acid 32449-15-3, N-1-Naphthalenyl-1,2-ethanediamine monohydrochloride  
RL: ARG (Analytical reagent use); DEV (Device component use); ANST (Analytical study); USES (Uses)  
(nitrogen dioxide detection by optical **sensor** with liquid-core-filled optical fiber waveguide containing indicator solution)  
IT 7732-18-5, Water, analysis  
RL: AMX (Analytical matrix); ANST (Analytical study)  
(organic compds. detection in water by optical **sensor** with liquid-core-filled optical fiber waveguide containing indicator solution)  
IT 7782-44-7, Oxygen, analysis  
RL: ANT (Analyte); ANST (Analytical study)  
(oxygen detection by optical **sensor** with liquid-core-filled optical fiber waveguide containing fluorescent indicator solution)  
IT 613-11-6, Leucomethylene blue 7758-89-6, Cuprous chloride 10049-05-5,

Chromous chloride  
 RL: ARG (Analytical reagent use); DEV (Device component use); ANST  
 (Analytical study); USES (Uses)  
 (oxygen detection by optical **sensor** with liquid-core-filled  
 optical fiber waveguide containing fluorescent indicator solution)  
 IT 10028-15-6, Ozone, analysis  
 RL: ANT (Analyte); ANST (Analytical study)  
 (ozone detection by optical **sensor** with liquid-core-filled  
 optical fiber waveguide containing indicator solution)  
 IT 81-88-9 13558-31-1 17372-87-1, Eosin  
 RL: ARG (Analytical reagent use); DEV (Device component use); ANST  
 (Analytical study); USES (Uses)  
 (ozone detection by optical **sensor** with liquid-core-filled  
 optical fiber waveguide containing indicator solution)  
 IT 7722-84-1, Hydrogen peroxide, analysis  
 RL: ANT (Analyte); ANST (Analytical study)  
 (hydrogen peroxide detection by optical **sensor** with  
 liquid-core-filled optical fiber waveguide containing indicator solution)  
 RN 7722-84-1 HCAPLUS  
 CN Hydrogen peroxide (H2O2) (9CI) (CA INDEX NAME)

HO-OH

RE.CNT 2 THERE ARE 2 CITED REFERENCES AVAILABLE FOR THIS RECORD  
 ALL CITATIONS AVAILABLE IN THE RE FORMAT

L45 ANSWER 8 OF 16 HCAPLUS COPYRIGHT 2004 ACS on STN  
 AN 1999:524601 HCAPLUS  
 DN 131:283407  
 TI Field method for monitoring blood glucose in beef cattle  
 AU Rumsey, T. S.; Kahl, S.; Elsasser, T. H.  
 CS Growth Biology Laboratory, Livestock and Poultry Sciences Institute,  
 Agricultural Research Service, USDA, Beltsville, MD, 20705-2350, USA  
 SO Journal of Animal Science (Savoy, Illinois) (1999), 77(8), 2194-2200  
 CODEN: JANSAG; ISSN: 0021-8812  
 PB American Society of Animal Science  
 DT Journal  
 LA English  
 AB The purpose of this study was to determine the applicability of the Accu-Chek  
 Easy (ACE) human self-monitoring system for monitoring glycemic status in  
 cattle. The ACE method was compared with the Yellow Springs Instrument  
 (YSI) anal. laboratory method in two studies. A preliminary study (62 samples)  
 and a primary study (434 samples) involved a nine-fold range and a 10-fold  
 range, resp., of glucose concns. obtained during the acute phase response  
 of growing beef cattle to injections of varying dosages of endotoxin. The  
 ACE monitoring method compared with the YSI anal. method resulted in  
 similar patterns of glucose concentration change, similar ranking of glucose  
 means across endotoxin dosages during hyper- and hypoglycemia, and a close  
 relationship between paired YSI and ACE concns. from common samples. The  
 ACE method identified all nine animals that displayed hypoglycemic  
 distress during the acute phase response to endotoxin injection. The  
 relationship between the YSI anal. method and the ACE monitoring method  
 was found to be nonlinear ( $YSI = -38.2 + 13.6 \cdot ACE \cdot 50$ ;  $R^2 = .99$ ;  
 $Sy \cdot x = 7.3 \text{ mg/dL}$ ), and the use of this equation to predict YSI  
 values from ACE values in an independent data set resulted in linearity  
 when YSI was regressed on the predicted YSI values ( $YSI =$

-.78+1.00·Predicted YSI; R2 =.87; Sy·x = 6.9 mg/dL). Even though variation seemed greater for ACE than for YSI, we concluded that a system developed for human self-monitoring of blood glucose, such as the ACE, can be used to monitor the glycemic status of cattle.

CC 9-2 (Biochemical Methods)  
Section cross-reference(s): 14

ST blood glucose beef cattle YSI ACR method field monitoring

IT Colorimetry  
(ACE method; field method for monitoring blood glucose in beef cattle)

IT Glucose **sensors**  
(YSI method; field method for monitoring blood glucose in beef cattle)

IT Lipopolysaccharides  
RL: ADV (Adverse effect, including toxicity); BIOL (Biological study)  
(endotoxin; field method for monitoring blood glucose in beef cattle)

IT Toxins  
RL: ADV (Adverse effect, including toxicity); BIOL (Biological study)  
(endotoxins, from E.coli; field method for monitoring blood glucose in beef cattle)

IT Biosensors  
(enzymic; field method for monitoring blood glucose in beef cattle)

IT Acute-phase response  
Blood analysis  
Cattle  
Hyperglycemia  
Hypoglycemia  
(field method for monitoring blood glucose in beef cattle)

IT 50-99-7, D-Glucose, analysis  
RL: ANT (Analyte); ANST (Analytical study)  
(blood; field method for monitoring blood glucose in beef cattle)

IT 50-99-7, D-Glucose, analysis  
RL: ANT (Analyte); ANST (Analytical study)  
(field method for monitoring blood glucose in beef cattle)

IT 7722-84-1, Hydrogen peroxide, uses 9001-37-0, Glucose oxidase  
13746-66-2, Potassium **ferricyanide** 15244-10-7, Ferric sulfate hydrate  
RL: ARG (Analytical reagent use); ANST (Analytical study); USES (Uses)  
(field method for monitoring blood glucose in beef cattle)

IT 7722-84-1, Hydrogen peroxide, uses  
RL: ARG (Analytical reagent use); ANST (Analytical study); USES (Uses)  
(field method for monitoring blood glucose in beef cattle)

RN 7722-84-1 HCAPLUS

CN Hydrogen peroxide (H2O2) (9CI) (CA INDEX NAME)

HO-OH

RE.CNT 6 THERE ARE 6 CITED REFERENCES AVAILABLE FOR THIS RECORD  
ALL CITATIONS AVAILABLE IN THE RE FORMAT

L45 ANSWER 9 OF 16 HCAPLUS COPYRIGHT 2004 ACS on STN  
AN 1999:262104 HCAPLUS  
DN 130:275937  
TI Gas **sensors** based on optical properties of liquid-core-filled optical fiber waveguides  
IN Dasgupta, Purnendu K.; Liu, Su Yi; Fein, Harry  
PA World Precision Instruments, Inc., USA  
SO Eur. Pat. Appl., 15 pp.

KATHLEEN FULLER EIC 1700 REMSEN 4B28 571/272-2505

CODEN: EPXXDW

DT Patent  
LA English  
FAN.CNT 2

|      | PATENT NO.   | KIND | DATE     | APPLICATION NO. | DATE     |
|------|--|------|----------|-----------------|----------|
| PI   | EP 909946  | A2   | 19990421 | EP 1998-308481  | 19981016 |
|      | EP 909946  | A3   | 19990811 |                 |          |
|      | EP 909946  | B1   | 20040623 |                 |          |
|      | R: AT, BE, CH, DE, DK, ES, FR, GB, GR, IT, LI, LU, NL, SE, MC, PT, IE, SI, LT, LV, FI, RO  |      |          |                 |          |
|      | US 6011882   | A    | 20000104 | US 1997-951254  | 19971016 |
|      | US 6016372   | A    | 20000118 | US 1998-55865   | 19980406 |
| PRAI | US 1997-951254   | A    | 19971016 |                 |          |
|      | US 1998-55865  | A    | 19980406 |                 |          |
| AB   | A gas-permeable liquid-impermeable optical fiber waveguide containing a liquid core is used as a probe for the detection or measurement of a chemical compound, in which the waveguide core is filled with a light-transmitting reagent that undergoes a change in optical characteristics when exposed to the chemical compound. The optical fiber waveguide wall has a refractive index of <1.33. The chemical specie, if present, will permeate through the waveguide wall and react with or be absorbed in the core liquid. The waveguide typically contains Teflon AF 2400 [4,5-difluoro-2,2-bis(trifluoromethyl)-1,3-dioxole-tetrafluoroethene copolymer] as the waveguide material. Some examples of types of gases that can be detected include CO <sub>2</sub> , O <sub>2</sub> , CO, H <sub>2</sub> S, NO <sub>2</sub> , NH <sub>3</sub> , O <sub>3</sub> , H <sub>2</sub> O <sub>2</sub> , Cl <sub>2</sub> , concentrated acids, and detection of organic compds. in water. |      |          |                 |          |
| IC   | ICM G01N021-05   |      |          |                 |          |
|      | ICS G01N021-77; G02B006-20   |      |          |                 |          |
| CC   | 79-2 (Inorganic Analytical Chemistry)  |      |          |                 |          |
| ST   | optical fiber waveguide gas <b>sensor</b> ; core optical fiber waveguide gas <b>sensor</b>   |      |          |                 |          |
| IT   | Optical gas <b>sensors</b><br>Optical gas <b>sensors</b><br>(fiber-optic; gas <b>sensor</b> based on optical properties of liquid-core-filled optical fiber waveguides)  |      |          |                 |          |
| IT   | Optical gas <b>sensors</b><br>(gas <b>sensor</b> based on optical properties of liquid-core-filled optical fiber waveguides)   |      |          |                 |          |
| IT   | Fiber optic <b>sensors</b><br>Fiber optic <b>sensors</b><br>(gas; gas <b>sensor</b> based on optical properties of liquid-core-filled optical fiber waveguides)  |      |          |                 |          |
| IT   | Absorption spectroscopy<br>Colorimetry<br>Luminescence spectroscopy<br>Raman spectroscopy<br>(in gas sensing; gas <b>sensor</b> based on optical properties of liquid-core-filled optical fiber waveguides)  |      |          |                 |          |
| IT   | Acids, analysis<br>RL: ANT (Analyte); ANST (Analytical study)<br>(inorg., concentrated, detection of; gas <b>sensor</b> based on optical properties of liquid-core-filled optical fiber waveguides)  |      |          |                 |          |
| IT   | Fluorometry<br>(oxygen-quenched, in gas sensing; gas <b>sensor</b> based on optical properties of liquid-core-filled optical fiber waveguides)   |      |          |                 |          |
| IT   | 108-95-2, Phenol, uses<br>RL: ARG (Analytical reagent use); RCT (Reactant); ANST (Analytical study);   |      |          |                 |          |

- RACT (Reactant or reagent); USES (Uses)  
 (ammonia indicator; gas **sensor** based on optical properties of  
 liquid-core-filled optical fiber waveguides)
- IT 10025-98-6, Potassium tetrachloropalladate(II) 222159-57-1, Silver  
 p-sulfoaminobenzoate  
 RL: ARG (Analytical reagent use); RCT (Reactant); ANST (Analytical study);  
 RACT (Reactant or reagent); USES (Uses)  
 (carbon monoxide indicator; gas **sensor** based on optical  
 properties of liquid-core-filled optical fiber waveguides)
- IT 119-93-7, o-Tolidine 34314-06-2, Tetramethylbenzidine  
 RL: ARG (Analytical reagent use); RCT (Reactant); ANST (Analytical study);  
 RACT (Reactant or reagent); USES (Uses)  
 (chlorine indicator; gas **sensor** based on optical properties  
 of liquid-core-filled optical fiber waveguides)
- IT 124-38-9, Carbon dioxide, analysis 630-08-0, Carbon monoxide, analysis  
 7664-41-7, Ammonia, analysis 7722-84-1, Hydrogen peroxide  
 (H2O2), analysis 7782-44-7, Oxygen, analysis 7782-50-5, Chlorine,  
 analysis 7783-06-4, Hydrogen sulfide, analysis 10028-15-6, Ozone,  
 analysis 10102-44-0, Nitrogen dioxide, analysis  
 RL: ANT (Analyte); ANST (Analytical study)  
 (detection of; gas **sensor** based on optical properties of  
 liquid-core-filled optical fiber waveguides)
- IT 521-31-3, Luminol 13746-66-2, Potassium **ferricyanide**  
 RL: ARG (Analytical reagent use); RCT (Reactant); ANST (Analytical study);  
 RACT (Reactant or reagent); USES (Uses)  
 (hydrogen peroxide indicator; gas **sensor** based on optical  
 properties of liquid-core-filled optical fiber waveguides)
- IT 14402-89-2, Sodium nitroprusside  
 RL: ARG (Analytical reagent use); RCT (Reactant); ANST (Analytical study);  
 RACT (Reactant or reagent); USES (Uses)  
 (hydrogen sulfide and ammonia indicator; gas **sensor** based on  
 optical properties of liquid-core-filled optical fiber waveguides)
- IT 63-74-1, Sulfanilamide 121-57-3 32449-15-3, 1,2-Ethanediamine,  
 N-1-naphthalenyl-, monohydrochloride  
 RL: ARG (Analytical reagent use); RCT (Reactant); ANST (Analytical study);  
 RACT (Reactant or reagent); USES (Uses)  
 (nitrogen dioxide indicator; gas **sensor** based on optical  
 properties of liquid-core-filled optical fiber waveguides)
- IT 613-11-6, Leucomethylene blue 7758-89-6, Cuprous chloride 10049-05-5,  
 Chromous chloride  
 RL: ARG (Analytical reagent use); RCT (Reactant); ANST (Analytical study);  
 RACT (Reactant or reagent); USES (Uses)  
 (oxygen indicator; gas **sensor** based on optical properties of  
 liquid-core-filled optical fiber waveguides)
- IT 81-88-9 13558-31-1 17372-87-1, Eosin  
 RL: ARG (Analytical reagent use); RCT (Reactant); ANST (Analytical study);  
 RACT (Reactant or reagent); USES (Uses)  
 (ozone indicator; gas **sensor** based on optical properties of  
 liquid-core-filled optical fiber waveguides)
- IT 76-59-5, Bromthymol blue 143-74-8, Phenol red  
 RL: ARG (Analytical reagent use); ANST (Analytical study); USES (Uses)  
 (pH indicator; gas **sensor** based on optical properties of  
 liquid-core-filled optical fiber waveguides)
- IT 37626-13-4, Teflon AF 2400  
 RL: DEV (Device component use); USES (Uses)  
 (waveguide; gas **sensor** based on optical properties of  
 liquid-core-filled optical fiber waveguides)
- IT 7722-84-1, Hydrogen peroxide (H2O2), analysis

RL: ANT (Analyte); ANST (Analytical study)  
(detection of; gas **sensor** based on optical properties of  
liquid-core-filled optical fiber waveguides)

RN 7722-84-1 HCAPLUS

CN Hydrogen peroxide (H2O2) (9CI) (CA INDEX NAME)

HO-OH

L45 ANSWER 10 OF 16 HCAPLUS COPYRIGHT 2004 ACS on STN

AN 1998:771500 HCAPLUS

DN 129:350266

TI A novel assembly for perfluorinated ion-exchange membrane-based  
**sensors** designed for electroanalytical measurements in  
nonconducting media

AU Toniolo, Rosanna; Comisso, Nicola; Bontempelli, Gino; Schiavon, Gilberto;  
Sitran, Stefano

CS Department Chemical Sciences Technology, University Udine, Udine, I-33100,  
Italy

SO Electroanalysis (1998), 10(14), 942-947  
CODEN: ELANEU; ISSN: 1040-0397

PB Wiley-VCH Verlag GmbH

DT Journal

LA English

AB A perfluorinated ion-exchange membrane-based **sensor** suitable for  
electroanal. measurements in electrolyte-free media is described, which  
was assembled following a novel design enabling an easier preparation  
procedure. It was fabricated by inserting the terminal portion of a  
working Pt wire electrode into a Nafion tubing of suitable diameter and  
welding the wire thus wrapped to the bottom of a cell body by an  
insulating epoxy resin. The remainder upper part of the working electrode  
was covered by a Teflon tubing to avoid the elec. contact with the  
internal electrolyte introduced into the cell body, which was equipped  
with a counter and a reference electrode. As a result of this configuration,  
the actual working-electrode surface is the wire circumference contacted  
by the polyelectrolyte material at the bottom of the assembly which is  
exposed to the sample. The performance of this **sensor** was  
tested by cyclic voltammetry, amperometric monitoring and flow injection  
anal. for the electroanal. of a series of prototype analytes either  
dissolved in electrolyte-free water (H2O2, hydroquinone,  
**ferricyanide**, I- and Br-) or present in N2 atmospheres  
(triethylamine and O2). Detection limits for these analytes were estimated  
( $\sigma = 3$ ), together with the corresponding ranges within which the  
responses display a linear dependence on the analyte concentration. The novel  
assembly is suitable only for the anal. in electrolyte-free liquid samples,  
while for the anal. of gaseous atmospheres, especially for flowing gases,  
ion-exchange membrane **sensors** prepared by the more usual procedure  
based on the use of working electrode materials embedded into a moist  
polyelectrolyte membrane should be preferred.

CC 79-2 (Inorganic Analytical Chemistry)

Section cross-reference(s): 61

ST Nafion platinum polymer electrode membrane **sensor**; solid polymer  
electrolyte membrane **sensor**; amperometric **sensor**

Nafion platinum polymer electrode; voltammetric **sensor** Nafion  
platinum polymer electrode

IT **Sensors**

(amperometric; perfluorinated ion-exchange membrane-based **sensors** designed for electroanal. measurement in nonconducting media)

IT Epoxy resins, uses  
 RL: DEV (Device component use); PRP (Properties); USES (Uses)  
 (in perfluorinated ion-exchange membrane-based **sensors** designed for electroanal. measurement in nonconducting media)

IT Flow injection analysis  
 (perfluorinated ion-exchange membrane-based **sensors** designed as electroanal. FIA-detector in nonconducting media)

IT Membrane electrodes  
 Polyelectrolytes  
 (perfluorinated ion-exchange membrane-based **sensors** designed for electroanal. measurement in nonconducting media)

IT Solid electrolytes  
 (polymer; perfluorinated ion-exchange membrane-based **sensors** designed for electroanal. measurement in nonconducting media)

IT Ionomers  
 RL: DEV (Device component use); PRP (Properties); USES (Uses)  
 (polyoxyalkylenes, fluorine- and sulfo-containing; perfluorinated ion-exchange membrane-based **sensors** designed for electroanal. measurement in nonconducting media)

IT **Sensors**  
 (voltammetric **sensors**; perfluorinated ion-exchange membrane-based **sensors** designed for electroanal. measurement in nonconducting media)

IT 7727-37-9, Nitrogen, analysis  
 RL: AMX (Analytical matrix); ANST (Analytical study)  
 (determination by perfluorinated ion-exchange membrane-based **sensors** designed for electroanal. measurement in)

IT 7732-18-5, Water, analysis  
 RL: AMX (Analytical matrix); ANST (Analytical study)  
 (determination by perfluorinated ion-exchange membrane-based **sensors** designed for electroanal. measurement in electrolyte-free)

IT 121-44-8, Triethylamine, analysis 123-31-9, Hydroquinone, analysis  
 7722-84-1, Hydrogen peroxide, analysis 7782-44-7, Oxygen, analysis  
 13408-62-3, **Ferricyanide** 20461-54-5, Iodide, analysis  
 24959-67-9, Bromide, analysis  
 RL: ANT (Analyte); ANST (Analytical study)  
 (determination by perfluorinated ion-exchange membrane-based **sensors** designed for electroanal. measurement in nonconducting media)

IT 7440-06-4, Platinum, uses  
 RL: DEV (Device component use); PRP (Properties); USES (Uses)  
 (perfluorinated ion-exchange membrane-based **sensors** designed for electroanal. measurement in nonconducting media)

IT 7722-84-1, Hydrogen peroxide, analysis  
 RL: ANT (Analyte); ANST (Analytical study)  
 (determination by perfluorinated ion-exchange membrane-based **sensors** designed for electroanal. measurement in nonconducting media)

RN 7722-84-1 HCAPLUS

CN Hydrogen peroxide (H2O2) (9CI) (CA INDEX NAME)

HO-OH

L45 ANSWER 11 OF 16 HCAPLUS COPYRIGHT 2004 ACS on STN

KATHLEEN FULLER EIC 1700 REMSEN 4B28 571/272-2505



AN 1996:656965 HCAPLUS  
 DN 125:296650  
 TI Electrochemical system for rapid detection of biochemical agents that catalyze a redox potential change  
 IN Song, Herking; Hafeman, Dean G.  
 PA Molecular Devices Corporation, USA  
 SO U.S., 42 pp.  
 CODEN: USXXAM  
 DT Patent  
 LA English  
 FAN.CNT 1

|      | PATENT NO.     | KIND | DATE     | APPLICATION NO. | DATE     |
|------|----------------|------|----------|-----------------|----------|
| PI   | US 5567302     | A    | 19961022 | US 1995-483249  | 19950607 |
| PRAI | US 1995-483249 |      | 19950607 |                 |          |

AB The present invention relates to a system for detecting, in a reliable, precise and highly sensitive manner, biochem. agents such as enzymes that catalyze a redox potential change. One electrode is used to measure redox potential changes in an aqueous electrolyte containing the biochem. agents. Another electrode is used to deliver a feedback current to the electrolyte in response to measured changes in electrolyte redox potential. The amount of feedback current or charge delivered through the electrode to the electrolyte is sufficient in magnitude to maintain a constant redox potential. Quantitation of the amount of feedback current or charge necessary to maintain the constant redox potential may then be used to determine the amount of biochem. agents present. Alternatively, the redox potential need not be kept constant, but instead may be allowed to reach a new steady-state. Thus, the current, or charge, conducted by a feedback electrode to maintain a new steady-state potential in the presence of an enzymic reaction may be used to quantitate the amount of enzymic activity present. The present invention provides precision in the quantitation results, high sensitivity in enzyme detection, and a wider dynamic range for quantitation of the biochem. agent. The invention is especially useful for the determination of enzyme labels used in immunoassays, e.g.,  $\beta$ -D-galactosidase, horseradish peroxidase, alkaline phosphatase, and glucose oxidase.

IC ICM G01N027-26  
 NCL 205777500  
 CC 9-1 (Biochemical Methods)  
 Section cross-reference(s): 7, 15, 72, 76  
 ST immunoassay enzyme detn redox potential change; LAPS electrode array coulometric feedback system; semiconductor electrode electrolyte redox potential detn; light addressable potentiometric **sensor** enzyme detn  
 IT Coulometers  
 Electrodes  
 Electroluminescent devices  
 Electrolytes  
 Semiconductor devices  
 pH  
 (electrochem. system for detection of biomols. and enzymes that catalyze redox potential changes)  
 IT Albumins, analysis  
 RL: ANT (Analyte); ANST (Analytical study)  
 (electrochem. system for detection of biomols. and enzymes that catalyze redox potential changes)  
 IT Enzymes

RL: ANT (Analyte); CAT (Catalyst use); RCT (Reactant); ANST (Analytical study); RACT (Reactant or reagent); USES (Uses)  
 (electrochem. system for detection of biomols. and enzymes that catalyze redox potential changes)

IT Potentiometers  
 (light-addressable potentiometric **sensor**; electrochem. system for detection of biomols. and enzymes that catalyze redox potential changes)

IT Immunoassay  
 (enzyme, electrochem. system for detection of biomols. and enzymes that catalyze redox potential changes)

IT **Sensors**  
 (optical, electrochem. system for detection of biomols. and enzymes that catalyze redox potential changes)

IT Electric potential  
 (redox, electrochem. system for detection of biomols. and enzymes that catalyze redox potential changes)

IT 9001-37-0, Glucose oxidase 9001-37-0D, Glucose oxidase, biotinylated 9001-78-9, Alkaline phosphatase 9001-78-9D, Alkaline phosphatase, streptavidin conjugates 9003-99-0, Peroxidase 9003-99-0D, Peroxidase, biotinylated 9013-20-1D, Streptavidin, alkaline phosphatase conjugates 9031-11-2 9031-11-2D, biotinylated  
 RL: ANT (Analyte); ARG (Analytical reagent use); ANST (Analytical study); USES (Uses)  
 (electrochem. system for detection of biomols. and enzymes that catalyze redox potential changes)

IT 7722-84-1, Hydrogen peroxide, reactions 13408-62-3, **Ferricyanide** 13408-63-4, Ferrocyanide 54827-17-7, Benzidine, 3,3',5,5'-Tetramethyl-  
 RL: ARG (Analytical reagent use); RCT (Reactant); ANST (Analytical study); RACT (Reactant or reagent); USES (Uses)  
 (electrochem. system for detection of biomols. and enzymes that catalyze redox potential changes)

IT 7439-88-5, Iridium, analysis 7440-06-4, Platinum, analysis 7440-21-3, Silicon, analysis 7440-44-0, Carbon, analysis 7440-57-5, Gold, analysis 7631-86-9, Silicon oxide, analysis 12033-89-5, Silicon nitride, analysis  
 RL: ARU (Analytical role, unclassified); DEV (Device component use); ANST (Analytical study); USES (Uses)  
 (electrochem. system for detection of biomols. and enzymes that catalyze redox potential changes)

IT 7722-84-1, Hydrogen peroxide, reactions  
 RL: ARG (Analytical reagent use); RCT (Reactant); ANST (Analytical study); RACT (Reactant or reagent); USES (Uses)  
 (electrochem. system for detection of biomols. and enzymes that catalyze redox potential changes)

RN 7722-84-1 HCAPLUS  
 CN Hydrogen peroxide (H2O2) (9CI) (CA INDEX NAME)

HO-OH

L45 ANSWER 12 OF 16 HCAPLUS COPYRIGHT 2004 ACS on STN  
 AN 1996:114133 HCAPLUS  
 DN 124:192567  
 TI Studies of consumed chemiluminescence-based **sensors**

KATHLEEN FULLER EIC 1700 REMSEN 4B28 571/272-2505

AU Lu, Jian-Zhong  
 CS Department of Chemistry, Nanjing University, Nanjing, 210093, Peop. Rep. China  
 SO Huaxue Xuebao (1996), 54(1), 71-6  
 CODEN: HHHPA4; ISSN: 0567-7351  
 PB Kexue  
 DT Journal  
 LA Chinese  
 AB Six types of consumed chemiluminescence **sensors** for ascorbic acid, cyanide, Mn<sup>2+</sup>, Co<sup>2+</sup> and H<sub>2</sub>O<sub>2</sub> were developed. It was based on the new approach that all of the reagents involved in the chemiluminescence reactions were immobilized electrostatically on Amberlyst A-27 or D151 ion-exchange resin. The analytes of interest can be sensed directly by the reaction with the chemiluminescence reagents, which were eluted by Na<sub>3</sub>PO<sub>4</sub> or NaCl from the immobilization column prior to the chemiluminescence reaction. Not only these **sensors** have a wide linear range high sensitivity and simplicity of instrumentation, but also the immobilization methods of the chemiluminescence reagents are simple. They were applied successfully to the detns. of analytes in various samples.  
 CC 79-2 (Inorganic Analytical Chemistry)  
 Section cross-reference(s): 80  
 ST consumed chemiluminescence based **sensor**; ascorbic acid consumed chemiluminescence based **sensor**; hydrogen peroxide consumed chemiluminescence based **sensor**; manganese cobalt consumed chemiluminescence based **sensor**; cyanide consumed chemiluminescence based **sensor**  
 IT **Sensors**  
     (consumed chemiluminescence-based **sensors** for ascorbic acid, cyanide, Mn<sup>2+</sup>, Co<sup>2+</sup>, and H<sub>2</sub>O<sub>2</sub>)  
 IT 50-81-7, Ascorbic acid, analysis 57-12-5, Cyanide, analysis 7722-84-1, Hydrogen peroxide, analysis  
 RL: ANT (Analyte); ANST (Analytical study)  
     (consumed chemiluminescence-based **sensors** for)  
 IT 7439-96-5, Manganese, analysis 7440-48-4, Cobalt, analysis  
 RL: ANT (Analyte); ANST (Analytical study)  
     (divalent; consumed chemiluminescence-based **sensors** for)  
 IT 7440-50-8, Copper, analysis  
 RL: ARU (Analytical role, unclassified); DEV (Device component use); ANST (Analytical study); USES (Uses)  
     (divalent; in consumed chemiluminescence-based **sensors** for cyanide and H<sub>2</sub>O<sub>2</sub>)  
 IT 13746-66-2, Potassium **ferricyanide**  
 RL: ARU (Analytical role, unclassified); DEV (Device component use); ANST (Analytical study); USES (Uses)  
     (in consumed chemiluminescence-based **sensors** for ascorbic acid)  
 IT 521-31-3, Luminol 9074-22-0, Amberlyst A 27 163293-51-4, D 151  
 RL: ARU (Analytical role, unclassified); DEV (Device component use); ANST (Analytical study); USES (Uses)  
     (in consumed chemiluminescence-based **sensors** for ascorbic acid, cyanide, Mn<sup>2+</sup>, Co<sup>2+</sup>, and H<sub>2</sub>O<sub>2</sub>)  
 IT 7790-21-8, Potassium periodate  
 RL: ARU (Analytical role, unclassified); DEV (Device component use); ANST (Analytical study); USES (Uses)  
     (in consumed chemiluminescence-based **sensors** for manganese and cobalt)  
 IT 7440-47-3, Chromium, analysis

RL: ARU (Analytical role, unclassified); DEV (Device component use); ANST (Analytical study); USES (Uses)  
(trivalent; in consumed chemiluminescence-based **sensors** for H2O2)

IT 7722-84-1, Hydrogen peroxide, analysis  
RL: ANT (Analyte); ANST (Analytical study)  
(consumed chemiluminescence-based **sensors** for)  
RN 7722-84-1 HCAPLUS  
CN Hydrogen peroxide (H2O2) (9CI) (CA INDEX NAME)

HO--OH

L45 ANSWER 13 OF 16 HCAPLUS COPYRIGHT 2004 ACS on STN  
AN 1995:486710 HCAPLUS  
DN 123:131327  
TI Development of magnetic and electrochemical properties from the encapsulation of molecules in sol-gel glasses  
AU Lan, E. H.; Dave, B.; Dunn, B.; Valentine, J. S.; Zink, J. I.  
CS Dep. Materials Sci. Eng., Univ. California, Los Angeles, CA, 90024, USA  
SO Materials Research Society Symposium Proceedings (1995), 371(Advances in Porous Materials), 267-76  
CODEN: MRSPDH; ISSN: 0272-9172  
PB Materials Research Society  
DT Journal  
LA English  
AB The flexible solution chemical of the sol-gel process was used to encapsulate a wide variety of organic mols. and biomols. in the pores of inorg. matrixes. This paper describes two new types of sol-gel materials in which the dopant mols. induce specific magnetic and electrochem. properties. The encapsulation of ferritin, an iron storage protein, produces an optically transparent, paramagnetic sol-gel material. The size of the protein ( $\approx 100$  Å) makes this dopant among the largest mols. yet encapsulated by the sol-gel method. The 2nd material incorporates Fe(CN)<sub>6</sub> and exhibits mediated electron transport in the sol-gel matrix. The addnl. encapsulation of an enzyme (peroxidase or alc. dehydrogenase) leads to electrochem. detection of specific analytes via catalytic reactions.  
CC 79-2 (Inorganic Analytical Chemistry)  
Section cross-reference(s): 66, 80  
ST mol encapsulation sol gel glass **sensor**; ferritin encapsulation sol gel glass **sensor**; ferricyanide encapsulation sol gel glass **sensor**; enzyme encapsulation sol gel glass **sensor**; magnetic property mol encapsulation sol silica; electrochem property mol encapsulation sol silica  
IT Encapsulation  
**Sensors**  
(encapsulation of mols. in sol-gel glasses for preparation of magnetic and electrochem. **sensors**)  
IT Ferritins  
RL: ARU (Analytical role, unclassified); DEV (Device component use); ANST (Analytical study); USES (Uses)  
(magnetic and electrochem. properties from the encapsulation of mols. in silica gel)  
IT Silica gel, analysis  
RL: ARU (Analytical role, unclassified); NUU (Other use, unclassified); RCT (Reactant); ANST (Analytical study); RACT (Reactant or reagent); USES

(Uses)  
 (magnetic and electrochem. properties from the encapsulation of mols. in silica gel)

IT 9003-99-0, Peroxidase 9031-72-5, Alcohol dehydrogenase 13408-62-3, **Ferricyanide**  
 RL: ARU (Analytical role, unclassified); DEV (Device component use); ANST (Analytical study); USES (Uses)  
 (encapsulation of mols. in sol-gel glasses for preparation of electrochem. **sensors**)

IT 64-17-5, Ethanol, analysis **7722-84-1**, Hydrogen peroxide, analysis  
 RL: ANT (Analyte); ANST (Analytical study)  
 (encapsulation of mols. in sol-gel glasses for preparation of electrochem. **sensors** for)

IT 125495-77-4, Trimethylorthosilicate  
 RL: ARU (Analytical role, unclassified); DEV (Device component use); ANST (Analytical study); USES (Uses)  
 (in preparation of silicate sol for encapsulation of mols. for preparation of electrochem. **sensors**)

IT **7722-84-1**, Hydrogen peroxide, analysis  
 RL: ANT (Analyte); ANST (Analytical study)  
 (encapsulation of mols. in sol-gel glasses for preparation of electrochem. **sensors** for)

RN 7722-84-1 HCAPLUS  
 CN Hydrogen peroxide (H2O2) (9CI) (CA INDEX NAME)

HO-OH

L45 ANSWER 14 OF 16 HCAPLUS COPYRIGHT 2004 ACS on STN  
 AN 1993:598958 HCAPLUS  
 DN 119:198958  
 TI Microbiosensors for acetylcholine and glucose  
 AU Karube, Isao; Yokoyama, Kenji; Tamiya, Eiichi  
 CS Res. Cent. Adv. Sci. Technol., Univ. Tokyo, Tokyo, 153, Japan  
 SO Biosensors & Bioelectronics (1993), 8(3-4), 219-28  
 CODEN: BBIOE4; ISSN: 0956-5663  
 DT Journal  
 LA English  
 AB Microbiosensors based on carbon and platinum fibers are described. Carbon fibers were used to construct microelectrodes of 7  $\mu$ m diameter. Electrochem. operations for pre-electrolysis and measuring were examined for the highly sensitive determination of hydrogen peroxide. A triangular potential (-2 to +2 V vs. Ag/AgCl) was applied before measuring each pair of double pulses (first pulse: 750 mV; second pulse: 1100 mV). The determination limit was 0.1  $\mu$ M of hydrogen peroxide. The reproducible determination of hydrogen peroxide is possible even in samples containing albumin protein. The separation of hydrogen peroxide from ascorbic acid is also possible because the oxidation potential of ascorbic acid is different from that of hydrogen peroxide. An acetylcholine microsensor was fabricated by immobilizing acetylcholine esterase and choline oxidase on the carbon fiber by entrapment with poly(vinyl alc.)-quaternized stilbazole (PVA-SbQ). This **sensor**

gave a linear calibration plot for the range 0.1-1.0 mM with a linear correlation coefficient of 0.9842. Glucose oxidase (GOD) and glucose dehydrogenase (GDH) immobilized cylindrical platinum microelectrodes were fabricated, and their characteristics were evaluated, resp., by using 1,4-benzoquinone (BQ) and **ferricyanide** as electron mediators. Each enzyme was immobilized by using PVA-SbQ on a cylindrical microelectrode of 2  $\mu$ m diameter. A linear range in the calibration curve of the GOD-based glucose microsensor was observed to be wider than that obtained using a disk electrode of 1 mm diameter. The mediated response of the 2  $\mu$ m glucose **sensor** was compared with the response resulting from hydrogen peroxide detection. This result showed that a higher response and a wider linear range were observed with highly concentrated mediator. A much higher response of the GDH immobilized 2  $\mu$ m microelectrode was obtained when not only **ferricyanide** but also diaphorase was employed to reoxidize the NADH produced by the enzyme reaction. of GDH. The GDH-based glucose microsensor was found to be unaffected by the concentration of dissolved oxygen.

- CC 9-7 (Biochemical Methods)
- Section cross-reference(s): 2
- ST acetylcholine detn micro biosensor; glucose detn micro biosensor;
- biosensor micro glucose acetylcholine
- IT Albumins, miscellaneous
- RL: MSC (Miscellaneous)
- (hydrogen peroxide determination with microbiosensor in relation to)
- IT Immobilization, biochemical
- (of enzymes, for acetylcholine and glucose determination with
- microbiosensor)
- IT Electrodes
- (bio-, enzyme, micro-, for acetylcholine and glucose determination)
- IT 50-99-7, Glucose, analysis 51-84-3, Acetylcholine, analysis
- RL: ANT (Analyte); ANST (Analytical study)
- (determination of, microbiosensor for)
- IT 7722-84-1, Hydrogen peroxide, analysis
- RL: ANT (Analyte); ANST (Analytical study)
- (determination of, with microbiosensor, acetylcholine and glucose
- determination in
- relation to)
- IT 13408-62-3, **Ferricyanide**
- RL: ANST (Analytical study)
- (glucose determination response with diaphorase-glucose dehydrogenase-
- immobilized microbiosensor enhancement by)
- IT 50-81-7, Ascorbic acid, miscellaneous
- RL: MSC (Miscellaneous)
- (hydrogen peroxide determination with microbiosensor in relation to)
- IT 9000-81-1, Acetylcholine esterase 9028-67-5, Choline oxidase
- RL: PROC (Process)
- (immobilization of, for acetylcholine determination with microbiosensor)
- IT 9001-37-0, Glucose oxidase 9028-53-9, Glucose dehydrogenase
- RL: PROC (Process)
- (immobilization of, for glucose determination with microbiosensor)
- IT 9079-67-8, Diaphorase
- RL: PROC (Process)
- (immobilization of, with glucose dehydrogenase for glucose determination
- with
- microbiosensor)
- IT 7722-84-1, Hydrogen peroxide, analysis
- RL: ANT (Analyte); ANST (Analytical study)
- (determination of, with microbiosensor, acetylcholine and glucose
- determination in

relation to)  
 RN 7722-84-1 HCAPLUS  
 CN Hydrogen peroxide (H2O2) (9CI) (CA INDEX NAME)

HO-OH

L45 ANSWER 15 OF 16 HCAPLUS COPYRIGHT 2004 ACS on STN  
 AN 1993:111787 HCAPLUS  
 DN 118:111787  
 TI Ion-metal and ion-selective electrode properties compared on the basis of the polyelectrode model  
 AU Ilyushchenko, M. A.; Mirkin, V. A.; Falkenstern, L. E.  
 CS Kazakh State Univ., Almaty, Kazakhstan  
 SO Sensors and Actuators, B: Chemical (1992), B10(1), 21-9  
 CODEN: SABCEB; ISSN: 0925-4005  
 DT Journal  
 LA English  
 AB The behaviors of ion-metal, ion-selective and film electrodes in solns. containing a redox system (Fe(CN)6<sup>3-</sup>/Fe(CN)6<sup>4-</sup>, Fe<sup>3+</sup>/Fe<sup>2+</sup>, quinone-hydroquinone, H2O2, ascorbic acid and dissolved O2) are compared. The electrodes used are made of Ag and Ag halides and chalcogenides. The film electrodes based on Ag chalcogenides and the metallic electrode had the same properties. The Ag-halide electrode is similar to the ion-selective ones. The results are interpreted from the point of view of electrochem. kinetics and mixed conductivity depending on the stoichiometric composition  
 CC 72-2 (Electrochemistry)  
 Section cross-reference(s): 79  
 ST electrode property redox system potentiometric **sensor**; silver halide chalcogenide electrode redox system; halide silver electrode redox system **sensor**; chalcogenide silver electrode redox system **sensor**; potential silver halide chalcogenide redox system; film electrode redox system potentiometric **sensor**; ion selective electrode redox system **sensor**  
 IT Silver chalcogenides  
 Silver halides  
 RL: PRP (Properties)  
 (electrodes, properties of, effect of redox systems on, potentiometric **sensors** in relation to)  
 IT Electrodes  
 (ion-metal, properties of, potentiometric **sensors** in relation to)  
 IT Electric potential  
 (of film electrodes and ion-selective electrode, response of **sensor** in relation to)  
 IT **Sensors**  
 (potentiometric, effect of redox systems on)  
 IT Redox agents  
 (properties of ion-metal and ion-selective electrodes in presence of, response of potentiometric **sensor** in relation to)  
 IT Electrodes  
 (film, properties of, potentiometric **sensors** in relation to)  
 IT Electrodes  
 (ion-selective, properties of, potentiometric **sensors** in relation to)

IT Electrodes  
 (potentiometric, properties of)

IT 7440-06-4, Platinum, properties 7440-22-4, Silver, properties  
 12002-99-2, Silver telluride (Ag<sub>2</sub>Te)  
 RL: PRP (Properties)  
 (elec. potential of film electrode of, in solution containing redox system,  
 response of potentiometric **sensor** in relation to)

IT 7785-23-1, Silver bromide (AgBr)  
 RL: PRP (Properties)  
 (elec. potential of film electrode of, with and without silver  
 telluride in solution containing redox system, response of potentiometric  
**sensor** in relation to)

IT 7783-90-6, Silver chloride (AgCl), properties  
 RL: PRP (Properties)  
 (elec. potential of film electrode of, with and without silver  
 telluride, in presence of hydrogen peroxide, response of potentiometric  
**sensor** in relation to)

IT 7783-96-2, Silver iodide (AgI)  
 RL: PRP (Properties)  
 (elec. potential of ion-selective electrode of, in presence of ascorbic  
 acid, response of potentiometric **sensor** in relation to)

IT 21548-73-2, Silver sulfide (Ag<sub>2</sub>S)  
 RL: PRP (Properties)  
 (elec. potential of ion-selective electrode with, in solution containing  
 redox  
 couple, response of potentiometric **sensor** in relation to)

IT 50-81-7, Ascorbic acid, properties 106-51-4, Quinone, properties  
 123-31-9, Hydroquinone, properties 7439-89-6, Iron, properties  
 7722-84-1, Hydrogen peroxide, properties 7782-44-7, Oxygen,  
 properties 13408-62-3, **Ferricyanide** 13408-63-4, Ferrocyanide  
 RL: PRP (Properties)  
 (properties of ion-metal and ion-selective electrodes in solution containing  
 redox system with, response of potentiometric **sensor** in  
 relation to)

IT 7722-84-1, Hydrogen peroxide, properties  
 RL: PRP (Properties)  
 (properties of ion-metal and ion-selective electrodes in solution containing  
 redox system with, response of potentiometric **sensor** in  
 relation to)

RN 7722-84-1 HCAPLUS  
 CN Hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>) (9CI) (CA INDEX NAME)

HO-OH

L45 ANSWER 16 OF 16 HCAPLUS COPYRIGHT 2004 ACS on STN  
 AN 1993:68892 HCAPLUS  
 DN 118:68892  
 TI Non-linear and pulse phenomena during hydrogen peroxide reduction at  
 chalcopyrite (photo)cathodes  
 AU Cattarin, S.; Tributsch, H.  
 CS IPELP, Padua, 35100, Italy  
 SO Electrochimica Acta (1993), 38(1), 115-22  
 CODEN: ELCAAV; ISSN: 0013-4686  
 DT Journal  
 LA English

KATHLEEN FULLER EIC 1700 REMSEN 4B28 571/272-2505



- AB The current-voltage curves of H<sub>2</sub>O<sub>2</sub> reduction in an alkaline medium recorded at CuFeS<sub>2</sub> and CuInSe<sub>2</sub> cathodes showed non-monotonic profiles, with a pronounced current wave and a region of neg. i/U slope. In the latter region, photocurrents inverted in sign are observed at CuInSe<sub>2</sub> electrodes. XP spectra taken on CuFeS<sub>2</sub> electrodes after polarization expts. show products of surface corrosion and (depending on emersion potential) changes in the oxidation state of Cu. On the basis of electrochem. and spectroscopic results, the current wave is attributed to activation of a catalytic mechanism of H<sub>2</sub>O<sub>2</sub> reduction involving Cu species. Current oscillations are observed when the polarization conditions are properly set. An "elec." anal. of the oscillatory phenomena is proposed, focused on the conditions of polarization control and resulting circuit (in)stability. At CuInSe<sub>2</sub>, illumination may be used as a key parameter to switch the oscillatory regime on and off or to trigger individual oscillations. The system may be considered to be a simple model device of a light **sensor** based on an electrode/electrolyte junction.
- CC 72-2 (Electrochemistry)  
Section cross-reference(s): 67, 74
- ST hydrogen peroxide redn electrochem photoelectrochem; chalcopyrite cathode photocathode hydrogen peroxide; oscillation current redn hydrogen peroxide; copper iron sulfide electrode; indium copper selenide electrode
- IT Photoconductivity and Photoconduction  
(of copper iron sulfide in presence of hydrogen peroxide)
- IT Reduction, electrochemical  
(of hydrogen peroxide on copper iron sulfide or copper indium selenide electrodes)
- IT Oscillating reaction  
(photoelectrochem. reduction of hydrogen peroxide on copper iron sulfide and copper indium selenide electrodes)
- IT Reduction, electrochemical  
(photochem., of hydrogen peroxide on copper iron sulfide or copper indium selenide electrodes)
- IT Reduction catalysts  
(photoelectrochem., copper species, for hydrogen peroxide)
- IT 7664-41-7, Ammonia, uses 7757-82-6, Disodium sulfate, uses 10043-35-3, Boric acid, uses  
RL: USES (Uses)  
(electrochem. reduction of hydrogen peroxide at copper iron sulfide or copper indium selenide in solution containing)
- IT 13408-62-3, **Ferricyanide**  
RL: PRP (Properties)  
(electrochem. reduction of hydrogen peroxide on copper iron sulfide electrode in presence of, oscillation in relation to)
- IT 12015-76-8, Copper iron sulfide (CuFeS<sub>2</sub>) 12018-95-0, Copper indium selenide CuInSe<sub>2</sub>  
RL: PRP (Properties)  
(electrochem. reduction of hydrogen peroxide on, oscillation in)
- IT **7722-84-1**, Hydrogen peroxide, reactions  
RL: RCT (Reactant); RACT (Reactant or reagent)  
(reduction of, electrochem., on copper iron sulfide and copper indium selenide electrodes)
- IT 1310-58-3, Potassium hydroxide, uses 7447-40-7, Potassium chloride, uses 7631-99-4, Sodium nitrate, uses  
RL: USES (Uses)  
(voltammetry of hydrogen peroxide on copper iron sulfide in solution containing)
- IT **7722-84-1**, Hydrogen peroxide, reactions  
RL: RCT (Reactant); RACT (Reactant or reagent)

(reduction of, electrochem., on copper iron sulfide and copper indium  
selenide electrodes)

RN 7722-84-1 HCAPLUS

CN Hydrogen peroxide (H2O2) (9CI) (CA INDEX NAME)

HO- OH

=>